

Pharmacy products

BD Tip Cap



BD cat no.	Description	Packaging
308341	Luer tip cap	10/tray 2,000/case
305819	BD Luer-Lok™ single sterile tip cap	200/box 2,000/case
305822	Luer tip cap bulk-sterile tray	50/tray 2,000/case

BD Cornwall™ disposable syringe system

BD cat no.	Description	Packaging
305224	10 mL fluid BD Cornwall™ dispensing syringe	10/case
305225	Syringe tip connector	50/case

BD AutoShield Duo™ Safety Pen Needle

The BD AutoShield Duo is the only safety pen needle available in the U.S. that passively covers both the front and back end needles—providing a high level of protection from contaminated needlesticks after use. The dual-ended passive safety design can provide a higher level of peace of mind and protection in a fast paced environment where errors are costly.



BD cat no.	Description	NDC no.	Packaging
329515	30 G x 3/16 in. (5 mm)	08290-3295-15	100/box



BD PhaSeal™ closed system drug transfer device

The BD PhaSeal system ensures leakproof drug transfer thanks to sealed, elastomeric, double membranes. The BD PhaSeal system meets the NIOSH and ISOPP definitions that a CSTD “mechanically prohibits the transfer of environmental contaminants into the system and the escape of hazardous drug or vapor concentrations outside the system.”



BD PhaSeal protectors

BD cat no.	Description	Fits to vial to neck size	Pressure equalization	Packaging
515100	Protector 14 (P14) drug vial access device	13 mm	20 mL	50/box 200/case
515102	Protector 21 (P21) drug vial access device	20 mm	20 mL	50/box 200/case
515103	Protector 21 (P21) multi-pack drug vial access device	20 mm	20 mL	5 each shelf pack 250/case
515104	Protector 28 (P28) drug vial access device	28 mm	60 mL	25/box 100/case
515105	Protector 50 (P50) drug vial access device	20 mm	50 mL	25/box 100/case
515106	Protector 50 (P50) multipack	20 mm	50 mL	5 each shelf pack 150/case
515117	Protector 55 (P55) drug vial access device	20 mm	50 mL	25 each shelf pack 100/case

BD PhaSeal injectors



BD cat no.	Description	Packaging
515003	Injector luer lock (N35) syringe safety device	50 each shelf pack 200/case
515005	Injector luer lock (N35) multipack, syringe safety device	5 each shelf pack 250/case

BD PhaSeal connectors



BD cat no.	Description	Packaging
515200	Connector luer lock (C35) IV line access device	50/box 200/case
515202	Connector luer lock (C45) IV line access device for use with needle free ports	50/box 200/case

BD PhaSeal administration set family



BD cat no.	Description	Packaging
515301	Secondary set, drip chamber (C60) IV bag access device	20/box 80/case
515304	Y-Site connector (C80) with extension tubing	30/box 120/case
515306	Infusion adapter (C100) IV bag access device	50/box 100/case
515307	Infusion adapter (C100) IV bag access device multipack	5 each shelf pack 250/case

BD PhaSeal accessories



BD cat no.	Description	Packaging
515400	Assembly fixture (M12)	1/case
515401	Syringe tray (M15)	80/case
515403	Infusion clamp (M25)	25/box 100/case
515404	Cap for injector membrane (M50)	50/box 200/case
515405	Hanger for IV-bag (M55)	40/box 160/case

**As of August 2012