

Single-use surgical products

Insufflation products and accessories

V. Mueller™ insufflation products include adapters, tubing and Veress needles, adaptable to all insufflator brands for tubing standardization. The rotating luer lock connector on the durable, non-collapsing tubing helps prevent twisting or disconnection from the insufflator. The tubing is available with bidirectional 0.1 or 0.3 micron hydrophobic filters to reduce the chance of contamination. It is sterile and not made with natural rubber latex.



BD cat no.	Description	Packaging
Single-use insufflation tubing		
LA101-RLL	0.1 micron filter, 10' length	20/box
LA103-RLL	0.3 micron filter, 10' length	20/box
Reusable insufflation tubing adapters		
LAUNV-LL	Snowden-Pencer™, Circon, Smith & Nephew and Olympus® insufflators	1/box
LASTK-LL	Stryker® insufflators	1/box
LALNV-LL	Linvatec® insufflators	1/box
LAWLF-LL	Richard Wolf™ insufflators	1/box
LASTZ-LL	KARL STORZ™ insufflators	1/box
Single-use Veress needle		
007.00	14 G, 120 mm	12/box
007.01	14 G, 150 mm	12/box

DeFogger™ anti-fog solution

Each bottle of DeFogger™ anti-fog solution comes with a nonabrasive, x-ray detectable foam pad that allows universal placement and will not damage the laparoscope. It is sterile and not made with natural rubber latex.



BD cat no.	Description	Packaging
LAF0G-FT	6 gm flip-top bottle	20/box

Dissectors and sponges

Round stick sponges

Round stick sponges provide excellent absorbency, and they can be used as prep sponges.



BD cat no.	Description	Packaging
23275-740	1 in. (2.5 cm), cotton	10/pack 25 packs/case
23275-745	1 ½ in. (3.8 cm), cotton	10/pack 25 packs/case

Flat stick sponges

Flat stick sponges provide excellent absorbency, and they can be used as prep sponges for surgery or with a FOERSTER clamp for deep sponge sticking.



BD cat no.	Description	Packaging
23275-750	1 ½ in. (3.8 cm), cotton	10/pack 25 packs/case

K-dissector sponges

K-dissector sponges can be used as longitudinal or folded sponges for dissection (they include a count holder for easy sponge counts).



BD cat no.	Description	Packaging
23275-410	¼ in. x ¾ in. (.6 cm x 1.4 cm), cotton twill	5/pack 100 packs/case

Round hard dissectors

Round hard dissectors, also known as beehives, can be used for blunt dissection. They can be grasped with forceps through the small hole in the center (they include a count holder for easy sponge counts).



BD cat no.	Description	Packaging
23275-440	¾ in. (.9 cm), cotton	5/pack 100 packs/case

Rosebud dissectors

Rosebud dissectors can be used for blunt dissection and are designed for use in delicate surgery (they include a count holder for easy sponge counts).



BD cat no.	Description	Packaging
23275-460	¼ in. (.6 cm), cotton	5/pack 100 packs/case

Peanut sponges

Peanut sponges can be used for soft tissue dissection and sponging. They can be used as is or folded (they include a count holder for easy sponge counts).



BD cat no.	Description	Packaging
23275-470	¾ in. (.9 cm), cotton	5/pack 100 packs/case

Cherry sponges

Cherry sponges provide excellent absorbency, and they can be used for general dissection (they include a count holder for easy sponge counts).



BD cat no.	Description	Packaging
23275-480	½ in. (1.3 cm), cotton	5/pack 100 packs/case

Cylindrical sponges

Cylindrical sponges are highly absorbent and designed for use in tonsil, adenoid and nasal surgeries.



BD cat no.	Description	Packaging
23276-250	½ in. x 1 ½ in. (1.3 cm x 3.8 cm), strung, cotton	5/pack 50 packs/case

**All dissectors and sponges are x-ray detectable