


Lighted Articulating Ear Curette™



Bionix Lighted Articulating Curette combines illumination, magnification and a 90° articulating tip to help healthcare providers take less time to see and safely remove ear canal passage obstructions.

► PROBLEM:

Cerumen plugs and foreign bodies in the ear are often difficult to remove. Practitioners lack the ability to visualize or work from behind the occlusion to remove it.

► SOLUTION:

The Lighted Articulating Ear Curette brilliantly illuminates the ear canal and provides an articulating tip that allows the practitioner to get behind the occlusion.

► RESULTS:

Direct visualization and an articulating ear curette eliminates the blind, top down sweeping motion that can abrade the ear canal improving safety and efficiency.



©2016 Bionix Development Corp. Patent Pending.

PRODUCT OFFERING

#2511 – Lighted Articulating VersaScoop® (3mm tip)

Includes 25 Articulating Ear Curette Tips,
1 Light Source and 1 Magnification Lens.

#2200 – Extra Light Source

#2265 – Pack of 5 Magnification Lenses

POTENTIAL REIMBURSEMENT (U.S. ONLY)

CPT Code / Description

69210 - For potential reimbursement of impacted cerumen.

69200 - For potential reimbursement of foreign body removal from the ear.

30300 - For potential reimbursement of foreign body removal from the nose.



Obelis s.a • Boulevard Général Wahis 53 • 1030 Brussels, BELGIUM
Tel: +(32) 2.732.59.54 • Fax: +(32) 2.732.60.03 • E-Mail: mail@obelis.net



Bionix Development Corp. • 5154 Enterprise Blvd. Toledo, OH, 43612, U.S.A.
Phone: 800-551-7096 • Fax: 800-455-5678 • www.Bionix.com

Phone: 800.551.6810
www.Bionix.com

BIONIX® GOING-BEYOND ►