


Lighted Placement Tool



Bionix Placement Tool combine illumination, magnification, and 5 depth guides to help Audiologists take less time to safely and accurately place a foam block deep in the ear.

► PROBLEM:


Placing a foam block deep in the ear canal is difficult without proper visualization.

► SOLUTION:

The Lighted Placement Tool combines the award winning technology of the Lighted Ear Curette™ with Magnification with a unique, bayonet design for an enhanced site line.

► RESULTS:

An elegantly slim, brilliantly lit and simply magnified device to promote proper visualization of the ear canal for more accurate foam block placement.



©2016 Bionix Development Corp. Patent Pending.

PRODUCT OFFERING

#2850 - Lighted Placement Tool

(1 Light Source, 1 Magnification Lens, 15 Lighted Placement Tips)

#2200 - Extra Light Source

#2265 - Pack of 5 Magnification Lenses

#2299 - Light Source Lanyard



Phone: 800.551.6810
www.Bionix.com

B BIONIX® GOING-BEYOND ►