

BLICKMAN



Tables

Adaptable work and prep surfaces

Includes:

- Instrument Tables
- Utility Tables
- Work Tables



OUR MISSION

Customer Focus

Exceed customer expectations through a commitment to continuous improvement of our products and services.

Quality Products

Offer the marketplace a broad range of high quality, cost effective, and durable products of superior craftsmanship.

Organizational Excellence

Strive for excellence in all we do by fostering a participative, team based environment, which nurtures the active involvement of all members of the organization.

Marketplace Recognition

Earn the respect, admiration, and loyalty of our people, customers, and suppliers by treating them with dignity and respect, embracing diversity, and developing and maintaining long term partnerships.

Collaborative Workplace

Encourage all employees to work safely and efficiently, contribute innovative ideas, and personally benefit from the ongoing success of the company.


BLICKMAN

THE BLICKMAN DIFFERENCE

Adjustable Tables

- Our adjustable tables are instrument tables designed to straddle the OR table using either manual or motorized adjustable height mechanisms to accommodate changes in the table position during neuro, cardiac, and orthopedic surgery, as well as to accommodate varying heights of personnel during general surgery.
- You can customize the tabletop dimensions to fit your differing surface area needs and the leg clearance dimensions to fit your differing OR table dimensions.
- Both models are made out of stainless steel and include casters for mobility and both models. The motorized table moves along telescoping aluminum lifting columns.

Key Product Features




Tables

Adjustable Tables


THE LIBERTY TABLE

Description: Our Liberty table was designed to straddle the OR table with adjustable height mechanisms. It travels effortlessly on antistatic swivel casters with brakes and utilizes a motorized height adjustment mechanism. Use a foot pedal or hand remote to operate the motorized adjustment. Two standard surface areas available.

Customize the surface area dimensions to fit your needs.



Standard Features



Foot Pedal



Hand Remote

MOTORIZED ADJUSTMENT:

Use the motorized foot operated/hand operated height adjustment mechanism.

Name	Liberty Table
Dimensions	41" W x 46" D (29" Leg Clearance)
Model #	7902SS-M
Item #	0157902000
Shipping Weight	315 lbs.
Dimensions	50 1/2" W x 24" D (38 3/8" Leg Clearance)
Model #	7901SS-M
Item #	0157901000
Shipping Weight	290 lbs.

Custom Dimensions Available Upon Request
Please email us at info@blickman.com


Tables

Adjustable Tables

THE GEORGETOWN TABLE

Description: Our Georgetown table was designed to straddle the OR table with a manual hand crank adjustable height mechanism. Travels easily on 3" swivel casters with brakes. Two standard surface areas available. Select extended arm design for electrophysiology procedures.

Customize the surface area dimensions to fit your needs.



Height adjustable from 36" to 56" off the ground.



Travels easily on four 3" swivel casters, two with brakes.

Standard Features



Hand Crank

MANUAL ADJUSTMENT:
Manually adjust the height of this table from either side with a detachable crank handle.

62" CRT table

This Georgetown table with the extended arm was designed for optimizing electrophysiology procedures.



Name	Georgetown Table
Dimensions	36" W x 20" D (28 5/8" Leg Clearance)
Model #	7892SS-0
Item #	0157892000
Shipping Weight	135 lbs.
Dimensions	48" W x 24" D (40 5/8" Leg Clearance)
Model #	7895SS-0
Item #	0157895000
Shipping Weight	160 lbs.
Dimensions	62" CRT Table (14" Overhang)
Model #	7895SS-O CRT
Item #	0157895062
Shipping Weight	190 lbs.

Custom Dimensions Available Upon Request
Please email us at info@blickman.com

Tables

Adjustable Tables

HYDRAULIC INSTRUMENT TABLE

Description: Our hydraulic instrument table is height adjustable to ensure compatibility with multiple OR table heights and healthcare personnel heights. Adjust its height with a manual hand crank from 30" - 42". It includes one lower shelf and stands stationary on glider feet tips. Select the surface area with dimensions that fit your needs.

Customize the surface area dimensions to fit your needs.



Name	Georgetown Table
Dimensions	48" x 24"
Model #	7835HYD
Item #	0177835000
Shipping Weight	145 lbs.
Dimensions	60" x 30"
Model #	7895SS-0
Item #	0177836030
Shipping Weight	188 lbs.

Custom Dimensions Available Upon Request
Please email us at info@blickman.com

Tables

Adjustable Tables

OVERBED TABLE

Description: Our overbed table comes in an attractive mark resistant laminate top in a walnut wood grain. It moves up or down with a simple one touch adjustment. Security features include a built-in column lock and four point base. The table travels smoothly on four 2" hooded casters.

Overbed Table Dimensions

Width: 30" x Height: 30" (max) x Depth: 15 1/2"



Name	Overbed Table
Model #	3400
Item #	2833400000
Shipping Weight	27 lbs.

Construction Specifications

Universal Dimensions	Adjustable heights to fit over OR table or hospital bed
Materials	304 stainless steel Liberty: Telescoping aluminum lifting columns. Overbed Table: Mar resistant laminate (tray), chrome plate steel (column & base)
Gauge	Liberty: 16GA Top; 11GA Sub Top; 14GA Base Georgetown: 18GA Top; 16GA Sub Top, 2" 16GA square tube legs
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	Liberty: 4" antistatic swivel casters with brakes Georgetown: 3" swivel casters, two with brakes
Notes	Liberty: Motorized adjustment mechanism Electrical Ratings: 100-120VAC, 15AMP, 60 Hz. Georgetown: Manual adjustment mechanism

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Custom Dimensions Available Upon Request
Please email us at info@blickman.com

THE BLICKMAN DIFFERENCE

Stationary Tables

- We offer a variety of stationary tables to fit any healthcare need. Our stationary tables include instrument tables, utility tables, and general work tables. Some models also come in MR conditional equivalents.
- Our instrument tables are reinforced with H-braces for stability or include an optional lower shelf for storing supplies. We offer standard dimensions or that can be customized to any size. We also have an optional nested set of up to six tables that can be ordered together or separately.
- Our utility tables come in a variety of drawer and shelf configurations to store and protect equipment and can include a workspace so you can store and prep your supplies on one efficient and transportable, setup.
- Our work tables are heavy duty with large surface areas designed for use in sterile and non-sterile preparation.

Key Product Features



Tables

Instrument Tables

NESTED TABLES

Description: Our nested instrument tables are space saving alternatives with a unique design that requires a mere 8 square feet of floor space, but yields 31 square feet of usable work space. All tables travel on 3" casters.



Name	Dimensions (WxHxD)	Model #	Item #	Ship Weight
Complete Set of 6	8 Sq. Ft. of Floor Space 31 Sq. Ft. of Work Space	7821SS	0137821000	255 lbs.

ORDER SEPARATELY

Largest	48" x 42" x 24"	7822SS	0137822000	56 lbs.
Second	44" x 40" x 22"	7823SS	0137823000	52 lbs.
Third	40" x 38" x 20"	7824SS	0137824000	51 lbs.
Fourth	36" x 36" x 18"	7825SS	0137825000	37 lbs.
Fifth	32" x 34" x 16"	7826SS	0137826000	34 lbs.
Smallest	28" x 32" x 14"	7827SS	0137827000	32 lbs.

MR options available, see MR section for more details



Tables

Instrument Tables

HOWARD & SAWYER TABLES

Description: Our instrument tables offer complete 18 gauge, 304 stainless steel construction. They are lightweight, yet extremely durable. Sanitary closed top help prevents bacterial build up. Sealed reinforcing channels provide rigidity. All tables are 34" tall, which is in the safety zone for reduced back stress, and travel on 3" casters. 16" clearance between shelves.

Sawyer Table



Howard Table



Table	H-Brace	Shelf
20" x 16"		
Model #	7840SS	7830SS
Item #	0127840000	0117830000
Shipping Weight	28 lbs.	32 lbs.
36" x 24"		
Model #	7843SS	7833SS
Item #	0127843000	0117833000
Shipping Weight	48 lbs.	56 lbs.
48" x 24"		
Model #	7845SS	7835SS
Item #	0127845000	0117835000
Shipping Weight	56 lbs.	69 lbs.
60" x 24"		
Model #	7846SS	7836SS
Item #	0127846000	0117836000
Shipping Weight	189 lbs.	200 lbs.

Options: Custom dimensions, three sided guardrails, and drawers (see Crescent table design).

MR options available, see MR section for more details



ANGULAR INSTRUMENT TABLE

Description: Choose our angular instrument table to have separate left and right work spaces. Add a lower shelf or drawers for additional storage. This 34" tall table is in the safety zone for reduced back stress, and it travels on 3" casters. Select the surface area with dimensions that fit your needs.



Name	Angular Instrument Table
Dimensions	60" x 15"
Model #	7862SSA
Item #	0147862000
Shipping Weight	120 lbs.
Dimensions	72" x 18"
Model #	7863SSA
Item #	0147863000
Shipping Weight	137 lbs.

Custom Dimensions Available Upon Request
Please email us at info@blickman.com

Tables

UTILITY TABLES

Description: Our utility tables are popular in the OR and throughout healthcare facilities as a general utility table. Our superior construction adds a supporting channel underneath the top shelf to enhance its strength. All of our utility tables are 34" tall, which is in the safety zone for reduced back stress, and travel easily on 3" casters.

Shelf Configuration Includes:

- ✓ Two shelves, 17" and 34" off the ground
- ✓ Four sided guardrails for protection

Drawer Configurations Include:

- ✓ A positive stop to prevent accidental drawer removal
- ✓ Ball bearing slides for easy sliding drawers
- ✓ Three sided guardrails for protection (optional four sided guardrails are available)

All of our utility tables are the same dimensions:

Width: 20" x Height: 34" x Depth: 16"



Lower Shelf Clearance:
16"

Lower Shelf Clearance:
11"

Lower Shelf Clearance:
5 1/2"

Lower Shelf Clearance:
4 1/4"

Lowest Drawer Clearance:
11 1/2"

# of Drawers	0	1	2	3	4
Model #	7852SS	7850SS	7754SS	7857SS	7859SS
Item #	0257852000	0217850000	0227754000	0237857000	0247859000
Shipping Weight	32 lbs.	46 lbs.	59 lbs.	72 lbs.	82 lbs.

Options: Feet: glides, rubber tips, or larger diameter casters, four sided guard rails on drawer models

Tables

Utility
Tables

CRESCENT TABLES

Description: Our crescent table has the same construction as the Howard and Sawyer Instrument tables with the addition of drawers and guard rails. This table includes a large work surface area (36"x 20"), two 4" deep drawers and one optional lower shelf 16" off the ground. This model is 34" tall, which is in the safety zone for reduced back stress, and travels easily on 3" casters.

Customize the surface area dimensions of these models to fit your needs.



Name	Crescent Table	
	H-Brace	Shelf
Model #	7856SS	7856SS-SH
Item #	0317856000	0317856003
Shipping Weight	78 lbs.	81 lbs.

Options: Feet: glides, rubber tips, or larger diameter casters

Tables

Work
Tables

WORK TABLES

Description: Our work tables are designed for use as prep tables. They are heavy duty with large surface areas designed for use in sterile and non-sterile environments. The smallest table dimensions uses four legs, and all other larger dimensions use six legs. All of our work tables are 36" tall, which is in the safety zone for reduced back stress, and come standard with bullet feet. Rubber tips, casters, and gliders are available as alternative options.

Customize the surface area dimensions of these models to fit your needs.



Surface Area	68" x 30"	80" x 36"	92" x 44"
		H-Brace	
Model #	CSDT6830	CSDT8036	CSDT9244
Item #	0336830000	0338036000	0339244000
Shipping Weight	210 lbs.	300 lbs.	390 lbs.
		Shelves	
Model #	CSDT6830S1	CSDT8036S1	CSDT9244S1
Item #	0336830100	0338036100	0339244100
Shipping Weight	260 lbs.	355 lbs.	435 lbs.

Options: Feet: glides, rubber tips, or larger diameter casters

Custom Designs Available Upon Request
Please email us at info@blickman.com



STATIONARY TABLE OPTIONS

Drawers

Description: For added storage, add drawers beneath any tabletop.



Drawer Options



Model #	Item #	Description	Shipping Weight
DH124S	2043424000	1 Drawer, 24" Wide Drawer Housing	35 lbs.
DH129S	2043429000	1 Drawer, 29" Wide Drawer Housing	50 lbs.
DH235S	2043435000	2 Drawers, 35" Wide Drawer Housing	55 lbs.
DH247S	2043447000	2 Drawers, 47" Wide Drawer Housing	60 lbs.

Custom Designs Available Upon Request
Please email us at info@blickman.com

STATIONARY TABLE OPTIONS

Guardrails

Description: For added placement security, add guardrails to any tabletop.



Feet Options

Rubber Tips | Casters | Gliders



Custom Designs Available Upon Request
Please email us at info@blickman.com

Tables

STATIONARY TABLE
Construction
Specifications

Construction Specifications

Universal Dimensions	Custom dimensions available on all tables
Materials	304 stainless steel
Gauge	Instrument Tables: 1 1/4" tube legs, 18GA table top Utility Tables: 1" Diameter 18GA tube legs, 18GA top & body Work Tables: 1 5/8", 16GA tube legs, 16GA table top & sub top
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	3" Casters
Notes	<ul style="list-style-type: none">MR conditional tables available where indicated

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Custom Dimensions/Designs Available Upon Request
Please email us at info@blickman.com

BLICKMAN

BLICKMAN

OUR PROMISE

For over 100 years, Blickman has been known throughout the healthcare industry for exceptional product quality. We pride ourselves on providing high-end and personalized service to ensure 100% customer satisfied. If you are not completely satisfied for any reason, we will make it right. All of our equipment is backed by our Limited Lifetime Warranty, so you can be certain that you have made the very best choice for your healthcare facility.

Call Blickman:
Toll Free: (800) 247-5070
Local: (973) 330-0557

www.blickman.com
info@blickman.com
500 U.S. Highway 46 East, Clifton, NJ 07011