

BLICKMAN



***Solution Stands/
Receptacles***

Roll, collect, dispose

Includes:

- Solution Stands
- Kickbucket
- Receptacle



OUR MISSION

Customer Focus

Exceed customer expectations through a commitment to continuous improvement of our products and services.

Quality Products

Offer the marketplace a broad range of high quality, cost effective, and durable products of superior craftsmanship.

Organizational Excellence

Strive for excellence in all we do by fostering a participative, team based environment, which nurtures the active involvement of all members of the organization.

Marketplace Recognition

Earn the respect, admiration, and loyalty of our people, customers, and suppliers by treating them with dignity and respect, embracing diversity, and developing and maintaining long term partnerships.

Collaborative Workplace

Encourage all employees to work safely and efficiently, contribute innovative ideas, and personally benefit from the ongoing success of the company.

BLICKMAN

THE BLICKMAN DIFFERENCE

- Our high quality, heavy-gauge, noise-reducing, solution stands and receptacles ensure your solutions and waste products are easily collected and disposed of in these easy-to-move, easy-to-clean product.
- Our solution stands and receptacles move effortlessly and quietly on 2" swivel casters, and are made out of our high-quality, easy-to-clean 304 stainless steel with MR conditional options available for both.
- Our solution stands are 33" high for compatibility with surgical table heights and are in the safe zone for reduced back stress, and our receptacles are low to the floor to pass under surgical tables and make it easy to drop waste into without reaching around the table.
- Our basins and pails are removable for both the solution stands and receptacles that are held in place by rubber bumpers that also help prevent noise from metal on metal contact.

Key Product Features



Solution Stands

Description: We offer two types of stainless steel solution stands: a single basin model and a double basin model. Our solution stands are 33" tall offering compatibility with surgical table heights and are in the safe zone for reduced back stress. Choose between a lower shelf or an H-brace for added storage and support. MR conditional models available.









Dimensions:

Width: 15" x Height: 33" x Depth: 15"



Dimensions:

Width: 29" x Height: 33" x Depth: 15"

| Name | Single Basin Solution Stand | | Double Basin Solution Stand | | |
|--|---|---|-----------------------------|---|------------|
| |  |  | Shelf |  | |
| | | | |  | |
|  | Model # | 7807SS | 7807MR | 7808SS | 7808MR |
| | Item # | 0717807000 | 0717807066 | 0727808000 | 0727808066 |
| | Shipping Weight | 20 lbs. | 20 lbs. | 33 lbs. | 33 lbs. |
| | | H-Brace | | | |
|  | Model # | 7807SS-HB | 7807MR-HB | 7808SS-HB | 7808MR-HB |
| | Item # | 0717807100 | 0717807166 | 0727808100 | 0727808166 |
| | Shipping Weight | 19 lbs. | 19 lbs. | 24 lbs. | 24 lbs. |

Solution Basin Also Sold Separately



| Name | Model # | Item # | Ship Weight |
|-----------------------------|---------|------------|-------------|
| 8.5 qt. Solution Basin Only | SB-1 | 9737370120 | 3 lbs. |

Receptacles

Description: Description: Our stainless steel kickbucket consists of a continuous encircling rubber bumper, which helps prevent tipping, wall scuff marks, and noise, a stainless steel cross brace, which provides excellent protection and a secure hold of a heavy gauge 13 qt. pail, and roller blade casters for mobility.

Parts Also Sold Separately



Kickbucket

13 qt. Pail



#KB-1



#7766NP

Parts Also Sold Separately



Receptacle

8.5 qt. Solution Basin





#SB-1



#8766NP

KICKBUCKET

| Name | Model # | Item # | Ship Weight |
|--|---------|------------|-------------|
| Lenox Kickbucket with Pail  | 7766SS | 0817766000 | 12 lbs. |
| Lenox Kickbucket with Pail  | 7766MR | 0817766066 | 12 lbs. |

ORDER SEPARATELY

| | | | |
|-----------------------------|--------|------------|--------|
| Lenox Kickbucket Frame Only | 7766NP | 0817766001 | 9 lbs. |
| 13 qt. Pail Only | KB-1 | 9701920000 | 4 lbs. |

RECEPTACLE

| Name | Model # | Item # | Ship Weight |
|------------------------------|---------|------------|-------------|
| Sponge Receptacle with Basin | 8766SS | 0828766000 | 11 lbs. |

ORDER SEPARATELY

| | | | |
|------------------------------|--------|------------|--------|
| Sponge Receptacle Frame Only | 8766NP | 0828766001 | 9 lbs. |
| 8.5 qt. Solution Basin Only | SB-1 | 9737370120 | 3 lbs. |



Solution Stands

Construction Specifications

Construction Specifications

| | |
|----------------------|---|
| Universal Dimensions | Solution Stands: 15", 29" x 33" x 15" Solution Basin: 8.5 qts. Pail: 13 qts. |
| Materials | 304 stainless steel |
| Gauge | Solution Stand Frame: 1" diameter tube, 16GA Solution Stand Shelf: 20GA Pail: 22GA; Basin: 14GA |
| Welding | All welding by heliarc method and exposed welds are ground smooth |
| Finish | Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish |
| Casters | 2" swivel casters |
| Notes | <ul style="list-style-type: none">MR conditional options available where indicated |

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

BLICKMAN

BLICKMAN

OUR PROMISE

For over 100 years, Blickman has been known throughout the healthcare industry for exceptional product quality. We pride ourselves on providing high-end and personalized service to ensure 100% customer satisfied. If you are not completely satisfied for any reason, we will make it right. All of our equipment is backed by our Limited Lifetime Warranty, so you can be certain that you have made the very best choice for your healthcare facility.

Call Blickman:
Toll Free: (800) 247-5070
Local: (973) 330-0557

www.blickman.com
info@blickman.com
500 U.S. Highway 46 East, Clifton, NJ 07011