

## OLSEN®

Olsen® – Pioneers in single-use bipolar forceps and electro-surgical instruments

# OPTIMAL CONTROL BETTER OUTCOMES

## DESIGNED FOR THERMAL CONTROL

### Olsen® Midas Touch™ Bipolar Forceps

The Midas Touch™ Non-Stick Coating maintains a smooth, highly conductive surface, enabling the use of lower power settings which may reduce tissue damage caused by thermal spread, eschar buildup and the production of noxious smoke.



- Minimizes tissue adherence – Reduces the need for repeated cleaning of the tips' surfaces
- Cermet coating is designed to maintain the instrument's non-stick properties

### Olsen® Bipolar Forceps

*Bipolar Forceps*



*Bayonet Forceps*



*Jewelers Forceps*



## DESIGNED FOR PRECISION

### Olsen® Micro Dissection Needles

*Ultra-Sharp Tungsten Wire*



## Why Olsen®?

- Market leader in bipolar forceps
- Option of proprietary Midas Touch™ non-stick coating
- Broad range of single-use bipolar forceps and micro dissection needles

# Olsen® Bipolar Forceps and Accessories

# Olsen® Micro Dissection Needles

OLSEN® MIDAS TOUCH BIPOLARS (5/box)	
20-1382IMT	Bayonet, Insulated, Smooth, Straight, 1.0 mm Tips, 191 mm (7 1/2 in)
20-1440IMT	Bayonet, Insulated, Smooth, Straight, 1.0 mm Tips, 210 mm (8-1/4 in)
20-1410IMT	Bayonet, Insulated, Smooth, Straight, 1.5 mm Tips, 210 mm (8- 1/4 in)
OLSEN® BIPOLARS (10/box)	
20-0421	Bayonet, Insulated, Irrigating, Smooth, Straight, 1.0 mm Tips, 210 mm (8- 1/4 in) (5/Bx)
20-0421R	Kit, Bayonet, Rotary, Insulated, Irrigating, Smooth, Straight, 1.0 mm Tips, 210 mm (8-1/4 in), (5/Bx)
20-0431	Bayonet, Insulated, Irrigating, Smooth, Straight, 1.5 mm Tips, 210 mm (8-1/4 in), (5/Bx)
20-1060*	Jewelers, Smooth, Straight Micro Tips, 114 mm (4-1/2 in)
20-1060I*	Jewelers, Insulated, Smooth, Straight Micro Tips, 114 mm (4-1/2 in)
20-1080*	Adson, Serrated, Straight, 1.0 mm Tips, 127 mm (5 in)
20-1080I*	Adson, Insulated, Serrated, Straight, 1.0 mm Tips, 127 mm (5 in)
20-1082I*	Adson, Insulated, Smooth, Straight, 1.0 mm Tips, 127 mm (5 in)
20-1090*	Semkin, Smooth, Straight, 0.75 mm Tips, 146 mm (5-3/4in)
20-1190KI	Kit, Semkin, , Insulated, Smooth, Straight, 1.0 mm Tips, 140 mm (5-1/2 in)
20-1360I	Bayonet, Insulated, Smooth, Straight, 0.75 mm Tips, 181 mm (7-1/8 in)
20-1370*	Bayonet, Smooth, Straight, 1.5 mm Tips, 181 mm (7-1/8 in)
20-1370I*	Bayonet, Insulated, Smooth, Straight, 1.5 mm Tips, 181 mm (7-1/8 in)
20-1380K	Kit, Bayonet, Smooth, Straight, 1.5 mm Tips, 191 mm (7-1/2 in)
20-1380I	Bayonet, Insulated, Smooth, Straight, 1.5 mm Tips, 191 mm (7-1/2 in)
20-1400KI	Kit, Bayonet, Insulated, Smooth, Straight, 0.75 mm Tips, 210 mm (8-1/4 in)
20-1410*	Bayonet, Smooth, Straight, 1.5 mm Tips, 210 mm (8-1/4 in)
20-1410I*	Bayonet, Insulated, Smooth, Straight, 1.5 mm Tips, 210 mm (8-1/4 in)
20-1415KI	Kit, Bayonet, Insulated, Smooth, Straight, 2.0 mm Tips, 210 mm (8- 1/4 in)
20-1440I*	Bayonet, Insulated, Smooth, Straight, 1.0 mm Tips, 210 mm (8-1/4 in)
20-1480I	Bayonet, Insulated, Smooth, Straight, 1.0 mm Tips, 267 mm (10-1/2 in)
ACCESSORIES & OTHER	
40-1102	12' (3.6m) Olsen Standard Twin Pin Bipolar Cord, Single Use, 10/Bx
40-2110	Monopolar Cord, 10 ft (3 m), Bovie Plug to Straight 4 mm Socket (Valleylab® Style), 10 ft, 3 m, Single Use, (10/bx)
40-3001R	12' (3.6m) Bipolar Irrigating Tubing and Bipolar Cord Set, For Rotary Pumps, Twin Pin, Single Use, 10/Bx

OLSEN® MICRO DISSECTION NEEDLES (5/box)	
97125BX	Needle Electrode, Ultra Sharp Tip, Insulated, Tungsten Wire, Bayonet Shaft, 15 cm (6 in)
97130BX	Needle Electrode, Tungsten Wire , Straight, Insulated, Ultra Sharp Tip, 3 cm (1-1/4 in)
97131BX	Needle Electrode, Straight, Ultra Sharp Tip, Insulated, Tungsten Wire, 4 cm (1-5/8 in)
97132BX	Needle Electrode, Ultra Sharp, Tungsten Wire , Insulated, Straight, 5 cm (2 in)
97133BX	Needle Electrode, , Tungsten Wire , Straight, Insulated, Ultra Sharp Tip, 6.4 cm (2- 1/2 in)
97134BX	Needle Electrode, Straight, Ultra Sharp Tip, Insulated, Tungsten Wire, 11 cm (4-5/8 in)
97135BX	Needle Electrode, Straight, Ultra Sharp Tip, Insulated, Tungsten Wire, 14 cm (5-1/2 in)
97136BX	Needle Electrode, Straight, Ultra Sharp Tip, Insulated, Tungsten Wire, 16.8 cm (6-5/8 in)
97141BX	Needle Electrode, 1 5/8" (4.1cm) , Tungsten Wire , 3.0mm Angled Tip, Insulated, Ultra Sharp Tip, 4 cm (1-5/8 in)
97143BX	Needle Electrode, 3.0mm Angled Tip, Ultra Sharp Tip, Insulated, Tungsten Wire, 5 cm (2 in)
97144BX	Needle Electrode, 3.0mm Angled Tip, Ultra Sharp Tip, Insulated, Tungsten Wire, 10 cm (4in)
97145BX	Needle Electrode, 3.0mm Angled Tip, Ultra Sharp Tip, Insulated, Tungsten Wire, 11.7 cm (4-5/8 in)
97147BX	Needle Electrode, 3.0mm Angled Tip, Ultra Sharp Tip, Insulated, Tungsten Wire, 16.8 cm (6-5/8 in)
97148BX	Needle Electrode, 3.0mm Angled Tip, Ultra Sharp Tip, Insulated, Tungsten Wire, 17.8 cm (7 in)
97150BX	Needle Electrode, 10mm Angled Shaft, Ultra Sharp Tip, Insulated, Tungsten Wire, 5 cm (2in)
97152BX	Needle Electrode, 10mm Angled Shaft, Ultra Sharp Tip, Insulated, Tungsten Wire, 12.4 cm (4-7/8 in)
97154BX	Needle Electrode, 10mm Angled Shaft, Ultra Sharp Tip, Insulated, Tungsten Wire, 17.5 cm (6 7/8 in)
97170BX	Needle Electrode, Curved Shaft, Ultra Sharp Tip, Insulated, Tungsten Wire, 10 cm (4 in)
97171BX	Needle Electrode, Curved Shaft, Ultra Sharp Tip, Insulated, Tungsten Wire, 15 cm (5-7/8 in)

\*Available as kit - Add K to PART CODE (Kit contains Forceps, Olsen® Standard Twin Pin Cord and Surgiwipe)