
Bovie Medical Corporation's electrosurgery business is now a part of



Suite of

OBGYN PRODUCTS

designed for your practice



Precision, versatility, and safety working seamlessly together

Expand patient care through more in-office procedures

Full suite of products offered in one convenient place







Bovie Better care for your patients from exam to procedure

The Next Generation of Clarity and Comfort in Colposcopy Exams

The Colpo-Master™ line offers superior, 3D, stereoscopic views with excellent clarity. Designed for maximum efficiency from exam to procedure, this suite of colposcopes features optimized ergonomics and a broad range of magnification.



Colpo-Master™ I Zoom Suspension Arm Colposcopes & Colpo-Master™ II Center-Post Mount Colposcopes

- 3D stereoscopic viewing
- · Zoom magnification range 4-27x with 20x eyepiece
- S-video and digital camera options available for connecting EMR platforms or monitor displays (trinocular head required)
- Dual light-intensity controls (green and white) for increased visibility of tissue variances
- Field of view ranges from 58.5mm at 4x to 8.5mm at 27x
- Binocular or camera-ready trinocular heads with 298mm working distance
- Independently adjust focus with dual diopter to eliminate eye strain

Colpo-Master™ III Triple Magnification Colposcopes

Capitalize on the versatility of a center-post-style design with 3-point magnification and semiautomatic focus.

- · 3D stereoscopic viewing
- Three magnification settings: 3.75x, 7.5x, and 15x
- Adjust magnification and maintain focus through knob rotation
- Flip-down green filter built into head of unit
- · Lightweight and easy to maneuver
- Dual diopter, independent, eye-focusing controls
- · 300mm working distance

Compare What's Best for Your Practice

	LED Illumination	Camera-ready Option	4-27x Zoom	Swing Arm	3-point Magnification	7-year Limited Warranty
Colpo-Master™ I	•	•	•	•		•
Colpo-Master™ II	•	•	•			•
Colpo-Master™ III	•				•	•

Procedure Safety and Reliability Your Practice Demands

Integrate the Bovie Specialist | PRO energy-based system with smoke evacuation to perform procedures like LLETZ/LEEP, wart removal, and other gynecologic procedures, allowing additional revenue streams by treating patients in-office.



Specialist | PRO High Energy-based Generator

Provide operating room-level technology to your patients in your office with this easy-to-use energy-based unit.


- Five modalities with up to 120 watts of power for a variety of tissue effects
- Monopolar and bipolar capabilities
- · Compatible with most industry-standard accessories
- Opportunity for additional revenue streams
- · 4-year warranty



Smoke Shark™ II Smoke Evacuator

Protect yourself, your staff, and patients from airborne pathogens and contaminants. Effectively removes hazardous smoke and particles from the operative site, reducing the risk of surgical site infection (SSI).

- 35-hour, long-life filter can reduce procedure costs by over 30%*
- Improves visibility of operative site
- Immediate evacuation of contaminated smoke particles
- Operates in conjunction with optional Bovie Remove Activation Switch for maximum filter life and surgical efficiencies
- Compatible with most energy-based generators



Disposable LLETZ/LEEP Loop

Achieve precise control during biopsies and other excisional procedures.

- 100% safety-tested insulation
- · Molded, ribbed safety grip
- · 21 sizes available



Disposable Ball Electrodes

Designed to control bleeding, fulgurate, and treat diseased tissue safely and efficiently.

- Safety grip insulators
- · Available in 5mm and 3mm balls
- Manufactured to the highest standards



Energy-based Mobile Stand

At your fingertips: keep surgical equipment organized and within reach at all times.

- Adjustable height extender
- · Instrument tray for instant access
- Bottom tray to store smoke evacuation equipment

Imagine working with more comfort, precision, confidence, and safety. Visit BovieMed.com for more information on Bovie's suite of women's health products.

Reimbursement Information

СРТ	Procedure	2015 National Fee¹
56605	Biopsy of vulva or perineum (separate procedure), one lesion	\$83.67
56820	Colposcopy of the vulva	\$114.06
56821	Colposcopy of the vulva with biopsy(s)	\$150.53
57420	Colposcopy of the entire vagina, with cervix, if present	\$120.49
57421	Colposcopy of the entire vagina, with biopsy(s)	\$160.54
57452	Colposcopy of the cervix including upper/adjacent vagina	\$110.84
57454	Colposcopy of the cervix with biopsy(s) of the cervix and endocervical curettage	\$155.53
57455	Colposcopy of the cervix with biopsy(s) of the cervix	\$145.52
57456	Colposcopy of the cervix with endocervical curettage	\$136.94
57460	Colposcopy of the cervix with loop electrode biopsy(s) of the cervix	\$286.04
57461	Colposcopy of the cervix with loop electrode conization of the cervix	\$323.94
57500	Biopsy, single or multiple, or local excision of lesion with or without fulguration (separate procedure)	\$129.79
57505	Endocervical curettage (not done as part of dilation and curettage)	\$103.33
58110	Endometrial sampling (biopsy) performed in conjunction with colposcopy (List separately in addition to code for primary procedure)	\$48.98
99170	Anogenital examination with colposcopic magnification in childhood for suspected trauma	\$161.25

Coding and reimbursement information is provided for informational purposes only. This content is not intended to instruct medical providers on how to bill for healthcare procedures or to describe proper coding, billing, and payment levels.

Reimbursement amounts ultimately paid to individual practitioners will differ due to geographic adjustments applied to national payment rates.

Beyond Energy-based Surgery

Bovie is a renowned industry leader offering a full suite of reliable and quality solutions for your practice to ensure positive patient outcomes while increasing profitability.

For product and purchasing information, please contact your representative at Sales@BovieMed.com

REFERENCE: 1. Wasserman Y, DMD, Medical Publishers, LTD. Physicians' Fee Reference. 2015;32. 169-172, 316.



Symmetry Surgical 3034 Owen Drive Antioch, TN 37013 USA