

## MI 300 AND MI 150

The best in high-value, space-saving, color-accurate LED exam lights. The MI 300 offers two brightness settings and its mobile base allows it to be moved anywhere easily. Both exam lights provide brilliant, color accurate illuminating and clarity. All lighting equipment (excluding normal replacement parts) is backed by a three (3) year limited warranty.



**XLDEM-150**



**XLDEM-300**

Dimensions:	
Unit Height	
MI 150	64.0" (1,625 mm)
MI 300	66.0" (1,676 mm)

Data:	
Light Source	LED
LED Life in Hours	50,000 (average)
Temperature	4,300° Kelvin
Lux MI 150	15,000 at 24" (60 cm) 25,000 at 18" (45.7 cm)
Lux MI 300	30,000 or 15,000 at 24" (60 cm) 50,000 or 25,000 at 18" (45.7 cm)
CRI	97
Operating Voltage	100 - 240 VAC 50/60 Hz
Light Head Input Voltage	24 VDC

MI-300 and MI-150	
XLDEM-150	MI 150 Exam Light Connectible Goose Neck w/Fixed Base (061223)
XLDEM-300	MI 300 Exam Light Flexible Goose Neck w/Rolling Base (061123)