

*From Dirty to Clean and Everything
In Between... It's BrandMax Clean!®*



TRI-CLEAN®

Tri-Clean by Brandmax Enzymatic Cleaners

Tri-Clean Enzymatic Cleaners provide superior cleaning power to effectively remove bio-burden. Our proprietary blend of enzymes and pleasant lavender fragrance will leave your instruments cleaner and work area refreshed! Combine Tri-Clean Ultrasonic Cleaners with Tri-Clean Enzymatic Cleaners for the *Ultimate Cleaning Solution!*



Tri-Clean Eco-System

PROTEASE
attacks proteins
& blood



AMYLASE
attacks starches,
sugars & carbohydrates

LIPASE
attacks body fats
& oils

Call Toll Free: 1-844-874-2532 or 404-418-8499

Email: CustomerService@BrandMax.us



Are Your Instruments Tri-Clean?

Take the *Tri-Clean* Challenge . . .

FEATURE	<i>Tri-Clean</i>	Competitor	BENEFIT
Proprietary blend of Protease, Amylase and Lipase enzymes	✓	?	Effective cleaning and removal of bio-burden
Proprietary surfactant blend to reduce surface tension	✓	?	Wetting action assists in lifting away bio-burden leaving your instruments cleaner
Neutral pH, non-corrosive and contains rust inhibitors	✓	?	Safe for all instruments, metals and surfaces
Soothing color and fragrance	✓	?	Clear visibility of soiled instruments and a refreshing aroma for your work environment
100% made in the USA	✓	?	From raw materials to finished product 100% made in the USA!

ITEM #	DESCRIPTION	DILUTION	QUANTITY
BEC32S	32 oz. Ready-To-Use Foam Spray	ready-to-use	12 quarts / case
BEC32PM	32oz Concentrate, Pre-measure	1 oz. per gallon	12 quarts / case
BEC1	Gallon, Concentrate	1 oz. per gallon	4 gallons / case

Want cleaner instruments?

Maintaining an elevated temperature of 110°F in your ultrasonic cleaner will increase enzyme activity and can improve infection protection!

Tri-Clean Eco-System

Heat enzymatic cleaner to 110°F

Powerful transducers to remove bio-burden



Triple enzymatic cleaner with Protease, Amylase and Lipase enzymes

Call Toll Free: 1-844-874-2532 or 404-418-8499

Email: CustomerService@BrandMax.us