



MILLENNIUM

SERIES Premium specialty seating with ergonomic design.



TASK

- 18.5" wide contour poly-foam 1 seat cushion; fully adjustable forward seat tilt
- Pneumatic height adjustment: 18" - 24.75"
- Adjustable backrest with synchronized 3 traversing movement of seat; black plastic shroud
- 4 Dual-wheel nylon casters
- 23" Five leg powder coated black aluminum base

Ergonomic design; seamless upholstery standard





#21435BA #21435B w/o arm rests

accessories - Millennium Series #03335 Qty 5 locking casters



- 18.5" wide contour poly-foam seat cushion
- Pneumatic height adjustment: 2 20.75" 31" w/glides 21.75" - 32" w/o glides
- Adjustable contoured backrest with black plastic shroud
- 4 18" foot ring
- 5 Dual-wheel nylon casters
- 6 23" Five leg powder coated black aluminum base

Ergonomic design; seamless upholstery standard



#21521B #21520B w/glides



SURGEON

- 18.5" wide contour poly-foam seat cushion
- Hands-free foot pneumatic height adjustment: 19" - 24"
- Adjustable contoured backrest with black plastic shroud
- Dual-wheel nylon casters; available with optional locking casters
- 23" Five leg powder coated black aluminum base

Ergonomic design; seamless upholstery standard






#21340

#21340V

#21340B #21340BV

backrestseamless/vacuum formed



The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company 13901 Main Street, Menomonee Falls, WI 53051 1.888.Brewer.1 www.brewercompany.com



CENTURY SERIES Preferred by medical facilities worldwide. Optional adjustable contoured backrest with black plastic shroud 16" wide, 4" thick poly-foam seat cushion Pneumatic or precision spinlift height adjustment Optional adjustable chrome foot ring Chrome hooded

VALUE PLUS Affordable without compromising style and reliability. SERIES

> 1 Optional adjustable contoured backrest with black plastic shroud

2 15" wide, 4" thick poly-foam seat cushion

3 Black plastic dual-wheel nylon casters

4 Pneumatic or spinlift height adjustment

5 Five-leg, 21" black ABS plastic composite base

17" – 22.75"

#22500V #22500B-C133

height adjustment

height range

#22500BV

#22500-C133 pneumatic

17" - 22.75"

#22400BV


#22400V

spinlift height adjustment

V - seamless/vacuum formed B - backrest C133 - Cal 133 upholstery


TRADITION

Economical exam stools that stand the test of time. SERIES



2 Sturdy four-leg, 19.5" chrome base

3 2" rubber wheel, ball bearing casters



Weight capcity for stools - 250 lbs.





Five-leg, 23" brushed aluminum

accent base

B
V - SEC
LCD - locking CC
G - glides
FR - foot ring
C133 - Cal 133 upholstery
UL - Ultraleather upholstery

backrest

height range

#11001BVFR #11001BFR-C133 #11001VFR

dual wheel casters

#11001B #11001BV #11001BLCD #11001BG* #11001B-C133 #1101B-UL

#11001V #11001LCD #11001G*

17.75" – 24.5"

#11001-C133 #11001-UL

pneumatic _ height adjustment