

Brewer
more possibilities

INSTRUMENT STANDS

INSTRUMENT STANDS

SERIES Ideal. Portable. Accessible.




Mayo Instrument Stand

- "California" style base for easy positioning
- 1" chrome plated steel tubing
- Two rubber wheel, ball bearing casters
- Removable stainless steel tray:
19.25" x 12.62" x .75"
- Base: 20" wheels, 19" diameter
- Height: 34" - 53"

#11100

accessories

#09412 replacement tray



Mobile Instrument Stand

- Adjustable height for comfortable use while sitting or standing
- 1" chrome plated steel tubing
- Stable five-leg cast aluminum base
- Removable stainless steel tray:
19.25" x 12.62" x .75"
- Dual wheel nylon casters
- Height: 29.25" - 48.75"

#43465

accessories

#09412 replacement tray



Electrosurgical Unit Stand

- Sturdy and stable for easy room-to-room portability and storage
- Adapts easily to most electrosurgical units
- Angled mounting panel for easy viewing
- Includes convenient cord and treadle hooks
- Stable five-leg cast aluminum base
- Dual wheel nylon casters
- Base: 22"
- Height: 37.75"

#43577

accessories

#43431 optional add-a-tray

#43433 optional clamp



Mayo Instrument Stand

- Designed for convenience during surgical procedures
- 1" chrome plated steel tubing
- U-shaped base minimizes interference with other equipment or patient
- Two rubber wheel, ball bearing casters
- Removable stainless steel tray: 19.25" x 12.62" x .75"
- Base: 19" wheels, 21" diameter
- Height: 32" - 50"

#43460

accessories

#09412 replacement tray

All Brewer instrument stands come with a 5-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company
13901 Main Street, Menomonee Falls, WI 53051
1.888.Brewer.1 www.brewercompany.com

Brewer
more possibilities