

Step Stools

Designed & built to support confident care.

Step Stool

Features rugged, steel construction with durable, black powder-coat finish and a skid-resistant rubber mat for safety. Dimensions: 9" H x 11.25" W x 14.25 L (22.8 cm x 28.6 cm x 36.2 cm). With handrail, total height is 35.25" (89.5 cm). Weight capacity: 350 lbs (158 kg).

Models

- 11200-1 - 1 Stool per box
- 11200-2 - 2 Stools per box
- 11220 - w/Hand rail

Accessories

- 99121 - Replacement crutch tips (set of 4)

Two-Step, Step Stool

Features heavy-duty, steel construction with durable, black powder-coat finish and a skid-resistant rubber mat for safety. Dimensions (1st step): 9" H x 12" W x 16 L (22.9 cm x 30.5 cm x 40.6 cm). Dimensions (2nd step): 7" H x 12" W x 16 L (17.8 cm x 30.5 cm x 40.6 cm). With handrail, total height is 44" (111.8 cm). Weight capacity: 600 lbs (272 kg).

Models

- 31200 - Step stool
- 31220 - Step stool w/hand rail

Accessories

- 50016-1 - Hand rail

Bariatric Step Stool

Features heavy-duty, steel construction with durable, black powder-coat finish and a skid-resistant rubber mat for safety. Dimensions: 9.25" H x 11.25" W x 14.25 L (23.5 cm x 28.6 cm x 36.2 cm). With handrail, total height is 34" (86.4 cm). Weight capacity: 600 lbs (272 kg).

Models

- 21200 - Step stool
- 21220 - Step stool w/hand rail

Accessories

- 99121 - Replacement crutch tips (set of 4)

Instrument Stands, Tables & Kick Bucket

Mobile Instrument Stand

Adjustable height from 29"-47.5" (73.6 cm-120.7 cm) for comfortable use sitting or standing. Removable stainless steel tray: 19.2" x 12.7" x 0.75" (48.8 cm x 32.3 cm x 1.9 cm).

Model 43465

Accessories

- 2102568 - Replacement tray

Mayo Instrument Stand

U-shaped base minimizes interference during procedures. Adjustable height from 31"-48.5" (78.7 cm-123.2 cm). Removable stainless steel tray: 19.2" x 12.7" x 0.75" (48.8 cm x 32.3 cm x 1.9 cm).

Model 43460

Accessories

- 2102568 - Replacement tray

Mayo Instrument Stand

"California" style base for easy positioning. Adjustable height from 34"-52" (86.4 cm-132.1 cm). Removable stainless steel tray: 19.2" x 12.7" x 0.75" (48.8 cm x 32.3 cm x 1.9 cm).

Model 11100

Accessories

- 2102568 - Replacement tray

Stainless Steel Mayo Stand

Features non-corrosive stainless steel construction. Adjustable height from 31"-48.5" (78.7 cm-123.2 cm). Removable stainless steel tray: 19.2" x 12.7" x 0.75" (48.8 cm x 32.3 cm x 1.9 cm).

Model 48461

Accessories

- 2102568 - Replacement tray

"U" Base Table

"U" base fits easily in restricted spaces. Offers convenient, spring-loaded height adjustment from 28.5"-43.25" (72.4 cm-109.9 cm). Top measures 15" x 30" (38.1 cm x 76.2 cm).

Models

- 11620 - w/Walnut grain laminate top & chrome base
- 11630 - w/Gray laminate top & gray epoxy finish on base

Kick Bucket w/Frame

Features non-corrosive stainless steel construction. Heavy-duty design includes 360° rubber bumper on frame 13 quart size (3.25 gallon).

Model 48711

- 48712 - Kick bucket only
- 48714 - Kick bucket frame only

Low "H" Base Table

Designed for low care beds. Features spring-loaded height adjustment. With casters: 19.375"-26.5" (49.2 cm-67.3 cm). With glides: 17.625"-24.5" (44.8 cm-62.2 cm). Top measures 15" x 30" (38.1 cm x 76.2 cm).

Model

- 11640 - w/Walnut grain laminate top & chrome base

"H" Base Table

Offers convenient, spring-loaded height adjustment from 28.5"-43.25" (72.4 cm-109.9 cm). Top measures 15" x 30" (38.1 cm x 76.2 cm).

Models

- 11600 - w/Gray laminate top & gray epoxy finish on base
- 11610 - w/Walnut grain laminate top & chrome base