

Additional Practice Seating

Ergonomic solutions for a wide range of needs.

Round out your clinic seating needs with rugged, reliable designs from Brewer. It's a choice your staff and your patients will appreciate.

Tradition Series

This economical exam stool is made to stand the test of time. The Tradition Series features spinlift height adjustment from 19.5"-27.5" (49.5 cm-69.9 cm), a sturdy four-leg chrome base, 14" (35.6 cm) wide, 3.5" (8.9 cm) thick poly-foam seat cushion and 2" (5.1 cm) rubber wheel, ball bearing casters.



Model 23051 in Arctic

Model 23051

Side Chair Series


Perfect for exam room, front office or waiting room environments, these durable chairs feature steel tube frames and kick-out rear leg design to help prevent wall dings and dents. Their contoured backrests and 2" (5.1 cm) thick poly-foam seat cushions offer plenty of comfort. Plus, they meet BIFMA (Business and Institutional Furniture manufacturer's Association) Standards.



Model 1200 in Ivy

Models

- 1200 - w/ sturdy plastic arm rests
- 1250 - w/o arm rests



Model 1250 Deep Sapphire upholstery

Blood Drawing Chair Series

Seat alignment positions patients securely on the chair, while innovative angled armrests and easy-to-use knob height adjustments provide optimal arm position for each patient's blood draw. Features include two standard add-on work surfaces (to accommodate prep materials, paperwork and vials). Chair also meets BIFMA Standards. An optional drawer (Model #101059) can be attached on either side of chair to keep prep materials nearby yet out of sight.



Model 1500 in Cabernet

Model 1500

Accessories

- 101059 - Drawer option

Upholstery Color Options

Standard upholstery color selections

For Brewer Side Chairs, Blood Drawing Chairs and Exam Stools



California Technical Bulletin 133 Upholstery Color Selections (CAL-133)

For Brewer Century Series Stools and Value Plus Series Pneumatic Stools

CAL TB 133 is a written legal standard outlining flammability test procedures for upholstered furniture. These vinyl color selections contain a flame blocking barrier that helps to prevent ignition, reflects heat and is flame retardant.



Standard upholstery color selections

For Brewer Side Chairs, Blood Drawing Chairs and Exam Stools

All Premium Plush Upholstery Colors are Ultrafabrics ULTRALEATHER™.



Colors may vary slightly depending on printing setting. For actual material swatches call 1•888•Brewer•1