




**Brewer**  
more possibilities

INFUSION PUMP STANDS

# INFUSION PUMP STANDS

**SERIES** Strength. Stability. Versatility.



## Large Infusion Pump Stand

- The most stable base available (26" diameter)
- Handles up to four infusion pumps and/or IV controllers (w/43430 and 43433 accessories)
- Tough six-leg epoxy coated steel base
- Chrome plated top post assembly
- 3" rubber wheel, ball bearing casters
- Exclusive "No Loss" knob
- Height: 52.25" - 102.75"
- Load capacity: 45 lbs.

#43408  
2 ram's horn hook

#43409  
4 ram's horn hook



## Short Wheel Infusion Pump Stand


- Heavy center weight affords stability while saving space (16" diameter)
- Tough five-leg epoxy coated steel base
- Chrome plated top post assembly
- 3" rubber wheel, ball bearing casters
- Exclusive "No Loss" knob
- Height: 54" - 90.5"
- Load capacity: 45 lbs.

#43416  
2 ram's horn hook

#43417  
4 ram's horn hook

## ACCESSORIES


The following accessories all attach with the #43433 clamp, which is sold separately



### Universal Clamp

- Fast, flexible attachment & positioning of accessories
- Set-screw locks accessory in place
- Secure grip; versatile adjustment


#43433



### Infusion Pump Frame

- Mount up to four pumps on a single stand
- Easy mounting, removal & adjustment of pump rack
- *Only suitable for use with #43408 and #43409 infusion pump stands*

#43430



### Steering Handle

- Allows easy & convenient mobility
- Non-slip black vinyl coating


#43428



### Oxygen Cylinder Holder

- Holds "D" and "E" size tanks
- Excellent for ambulatory patients

#43429



### Add-A-Tray

- Holds instruments or patient items
- Removable stainless steel tray
- Tray size: 19.25" x 12.62" x .75"

#43431

**All Brewer infusion pump products come with a 5-year warranty.**

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

**Call your local Brewer sales representative to learn about our complete product line.**

The Brewer Company  
13901 Main Street, Menomonee Falls, WI 53051  
1.888.Brewer.1 www.brewercompany.com

**Brewer**  
more possibilities