

CA100

SORRENTO

7 Section Therapeutic Table (Electric)



- 1 Year on Gas Struts, Castors, Upholstery
- 2 Years on Electronic Components
- 3 Years on Frame & Pillar



Prussian Blue Licorice Black Dim Grey Merlot (Red Wine) Creamy Beige



1 Motor



Hand Control



Foot Control



2 / 3 / 4 Sides Elevating Bar

Maximum User Weight
225 kgs / 496 lbs



Product Features

The CA100 – 7 Section Therapeutic Table has been our flagship since 2012. What made this model popular is its 15° Center Arching (Locked), 60° Reversible Head Section, Spring-Loaded Arm Rest with 120° Swivel and 87° Leg Incline capabilities all-in-one. Even better, the center joint can be UNLOCKED from the upper frame to make the leg section free-floating. A ramp & plateau positioning can be achieved as well as many other possible positions. It's simply versatile! Both arm rests can be opened and closed up to 120° as well as lifted and dropped. The two side wings just above the shoulder can be folded down when the user is lying down on her stomach. She can then extend her arms forward and rest on the arm rests comfortably.

- 60° Reversible Head Section (Length: 39cm / 15.35" inch)
- Head Section Incline: 25° / Recline: 35°
- 15° Center Arching (Locked) / 30° Center Arching (Unlocked)
- 87° Leg Section Incline (Length: 84cm / 33" inch)
- Spring-Loaded Arm Rest with 120° Swivel
- 1 Motorized Actuator for Hi/Lo Adj.
- 7 * Sections Upholstered Table Top
- Twin "Oversized" Pillars for Added Stability
- Individual "Leveling Foot" on All 4 Corners
- Retractable Castors at Both Ends
- Twin Gas Struts on All Sections for Solid Support

Technical Specifications

