

CA130

CHIROMA

8 Section Chiropractic Table (Electric)



- 1 Year on Gas Struts, Castors, Upholstery
- 2 Years on Electronic Components
- 3 Years on Frame & Pillar



Prussian Blue Licorice Black Dim Grey Merlot (Red Wine) Creamy Beige



3 Motors



Foot Control



2 Sides Elevating Bar

Maximum User Weight
225 kgs / 496 lbs

Product Features

This is a CA130 - 8 Section Chiropractic Table (Electric) with 4 Sections of Adjustable Tension-Drop at the (1) Head (2) Thorax and (3) Abdomen Sections.

As complex as it looks, the CHIROMA provides professional chiropractor with the simplicity of use. Since the main emphasis of chiropractic is on manual and manipulative therapy for the treatment of joint dysfunctions, especially those of the spinal column, CHIROMA serves as a "Force Extension" as well as "Force Absorbent" platform for chiropractors to perform various tentative tasks without the risks of collateral damage to the patient.

The CA130 Chiroma has the following features:

- Head Section is equipped with 2 motors for Hight and Tilting Adjustment
- 1* Joystick Control for the Head Section's Height and Tilting Adjustment
- 4* Sections are equipped with Adjustable Tension-Drop
- 2* Adjustable Arm Rests for absolute comfort & reliability
- 1* Motor for Height Adjustment
- Twin "Oversized" Pillars for Added Stability
- 2* Individual Leveling Foot at Head
- 2* Hauling Wheels at the Rear – Designed to be lifted and moved around by lifting just the Arm Rest



The adjustable tension drop can be initiated at all sections of our body:

- Cervical Drop
- Thoracic Drop
- Lumbar Drop
- Pelvic Drop

Technical Specifications

