


CA20^{DUO}

2 Section Therapeutic Table (Electric)



- 1 Year on Gas Struts, Castors, Upholstery
- 2 Years on Electronic Components
- 3 Years on Frame & Pillar



Prussian Blue Licorice Black Dim Grey Merlot (Red Wine) Creamy Beige



1 Motor



Hand Control



Foot Control



2/3/4 Sides Elevating Bar

Maximum User Weight
225 kgs / 496 lbs



Product Features

This is a CA20 - 2 Section Therapeutic Table with 1 Motor for Hi / Lo Adjustment. The 2 section therapeutic table has always been a popular model as the body & leg section is made into a single piece and provide sturdy support at all times.

- Hydraulic version (CA20H) is also available
- Equipped with 1 Motorized Actuator for Hi/Lo Adj.
- Twin "Oversized" Pillars for Added Stability
- Individual "Leveling Foot" on All 4 Corners
- Retractable Castors at Both Ends
- Twin Gas Struts for Solid Head Section Support

Technical Specifications

