


CA65 ^{TRIO}

3 Section Therapeutic Table (Electric)



- 1 Year on Gas Struts, Castors, Upholstery
- 2 Years on Electronic Components
- 3 Years on Frame & Pillar



Prussian Blue Licorice Black Dim Grey Merlot (Red Wine) Creamy Beige



1 Motor



Hand Control



Foot Control



2 / 3 / 4 Sides Elevating Bar

Maximum User Weight
225 kgs / 496 lbs



Product Features

The 3 Section Therapeutic Table is one of the most popular single product in existing product range. This is all due to its versatility in terms of adjustable head and foot drop sections.

- Hydraulic version (CA65H) is also available
- Equipped with 1 Motorized Actuator for Hi/Lo Adj.
- Twin "Oversized" Pillars for Added Stability
- Individual "Leveling Foot" on All 4 Corners
- Retractable Castors at Both Ends
- Twin Gas Struts for Solid Head and Leg Sections Support

Technical Specifications

