

CAPSAHEALTHCARE

Avalo® ACS Medication Cart Avalo® ACM Medication Cart

Smart, Secure & Simple

Capsa Healthcare brings its superior design expertise to the hospital market with the Avalo AC Medication Cart. This innovative product line features the latest design and serves as a platform for existing and emerging medication control technology.

Capsa has redefined medication cart security with a broad choice of lock systems including keyless entry with automatic relocking, mag stripe card access, keyless narcotic storage and cart user audit.

Best of all, AC Medication Carts offer huge capacity and intelligent solutions for medication management and organization.

AC Medication Cart Highlights & Options

- A Internal pass-thru waste receptacle
- **B** Optional expansive slide-out writing surface
- **C** Magnetic stripe, proximity or bar code security for cart access
- $\ensuremath{\mathsf{D}}$ Optional keyless access with auto-relocking
- **E** Integrated and flexible drawer divider system
- **F** Configurable medication storage (2 bin sizes & 3 drawer sizes)
- **G** Full selection of convenience accessories & options
- **H** Stable, trim-line footprint for simple maneuverability
- I Choice of sealed lead acid or LiFe battery system

Cassette & Bin Options

Designed with flexibility in mind, Capsa Healthcare's cassettes and bins give you a choice of 2 or 3 tier durable cassettes with the option of 5.5" or 8" bins.

2-Tier & 3-Tier Cassettes



The AC Medication Cart accommodates 2 & 3 tier cassettes to customize the storage configuration that best suits your requirements.

8" & 5.5" Bins



AC cassettes give you an option of 8" and 5.5" bins that offer a flexible divider system. Large bin lens permits protected and clear patient labeling.



Security Options

Capsa Healthcare has fitted AC Medication Carts with its proven and innovative automatic relocking system and the latest integrated medication cart security.

Lock System Choices

AutoLock™

Automatically secures a motionless unlocked medication cart with a userdefined relock time delay from 1-99 minutes. AutoLock permits up to 5,000 individual user codes to further enhance cart access tracking and security. *See available security upgrades.*

Keyless Access with Manual Relocking

Simply enter a 4-digit access code, push down the lock lever and the cart is ready for access. To secure the medication cart, simply push in the lock lever. Clear visual indication of lock status is apparent from a distance.

Key Access with Manual Relocking

Featuring the Best® Core Removable Lock System, this basic yet proven lock system provides your medication cart advanced security and the flexibility to change keys simply and easily.







To further optimize security, add a secondary access requirement such as ID badge card, bar code or proximity reader.

AutoLock[™] Security Upgrades

Keyless Narcotics Storage



The optional keyless narcotics storage feature greatly augments security to further optimize control and accountability.

Secondary Access



Select from a bar code, mag stripe or proximity card reader to serve as the primary or secondary access requirement

Audit Tracking



Authorized users can track cart activity, including user access and relocking for over 2,000 events.

Accessories & Options

Cart Light

Customize your configuration

with AC Medication Carts. Streamlined clinical workflow is assured with a full complement of convenience accessories and organizing supplies to facilitate patient care.



Clear Top Mat



IV Pole



Narcotics Storage Waste Container



Drawer Dividers





Storage Module



Sharps / Gloves



Mobile Computing

Specifications & Models

AC Medication Cart Features & Specifications

The Avalo AC Medication Cart line combines the latest medication cart design with plenty of capacity and simple maneuverability. Choose the drawer configuration and cassette size that best meets your requirements for highly-secure and organized medication control.

This solution is perfect for point-of-care medication administration and as a delivery workstation to augment your dispensing cabinet.

AC Cart Models

Avalo AC is available in 2 sizes and 2 heights:



Specifications, availability & components are subject to change without notice.

CE

Features

- Power status display on LCD
- Large retractable writing surface (optional)
- Internal pass-thru waste receptacle
- Modular bin / storage drawers
- Expansive storage drawers
- Slot-Lock™ security system secures all partially closed drawers and bins left open up to 1 inch
- Keyless access and auto-relocking with AutoLock option
- Wireless remote security control (optional)
- 5" double ball bearing casters (12.7cm)

Benefits

- Improved workflow and accuracy
- Safe and ergonomic design
- Advanced medication security and control
- Lightweight design for easy maneuverability
- Customized medication / supplies storage

Standard Cart Specifications

Base Dimensions:	ACS: ACM:	2.0.7.20 11	61cm x 66cm 61cm x 96.5cm
Work Surface Area:	ACS: ACM:		40.5cm x 46cm 40.5cm x 71cm
Work Surface Height:	Compa Standa	rd: 39.5" rd: 43"	100cm 109cm
Construction:	Formed steel, aluminum, heavy-duty polymer work surface / end panels		
Finish:	High-density polymer		
Wheels:	Swivel casters with 1 locking; tracking caster optional		

www.capsahealthcare.com | 800.437.6633

