

Avalo® Series Medical Cart

Emergency Cart-Convenient Packages

The Avalo Series Emergency Cart offers a standard of organization and function, with a host of features that provide reliable performance and promote improved workflow in emergency departments.

The Avalo Emergency Cart is available in blue or red in three different heights. Choose a base cart in Standard, Intermediate or Compact height, then add a convenient accessory package to meet the workflow and storage requirements for your facility.



Blue Emergency Cart
Standard Height
Shown with Accessory Package 2

Red Emergency Cart
Intermediate Height
Shown with Accessory Package 3



Quick-Access Lock

Avalo's unique breakaway locking handle offers optimal security and the simple access essential to timely emergency response. A swift downward swipe breaks the tamper evident seal, allowing access to all drawers.


Emergency Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement
- Quick-access lock offers secure, plus allows quick access for code response
- Smooth panels, work surface and seamless drawers promote simple cleaning for infection prevention
- Selection of accessories and organization systems
- 10-year limited warranty

Step 1

Choose a Base Cart Option:

Choose a cart, then add an accessory package to configure the cart for your specific requirements.



Cart Specifications

	Surface Height	Depth	Width
Standard	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
Intermediate	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
Compact	36" / 91 cm	24" / 61 cm	31" / 79 cm

Drawer Specifications

	Height	Depth	Width
Three Inch	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
Six Inch	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
Ten Inch	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

Specifications, availability & components are subject to change without notice.

Step 2

Choose an Accessory Package:

Choose your accessory package on the **next page**.

Step 2

Convenient Accessory Packages

Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

Avalo® Emergency Cart Accessory Package 1

AM-EM-ACCPK1



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole

Avalo® Emergency Cart Accessory Package 2

AM-EM-ACCPK2



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole



Defibrillator Shelf



3" Drawer Dividers (Set of 2)

Avalo® Emergency Cart Accessory Package 3

AM-EM-ACCPK3



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole



Defibrillator Shelf



3" Drawer Dividers (Set of 2)



Writing Surface



Sharps Container



Single Glove Box

Can't find what you're looking for?

If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.