



Avalo® Series Medical Cart

Procedure/Treatment Cart-Convenient Packages

The Avalo Series Procedure/Treatment Cart offers ample storage capacity to accommodate the largest bulk items to the smallest medical supplies.

The cart is adaptable to many different healthcare applications, including trauma, airway, dressing, cast, or general supply cart. A selection of proven lock systems provides secure and safe storage.



Procedure/Treatment Cart
Standard Height Shown with
Accessory Package 2

Procedure Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement and steering
- Smooth panels, work surface and seamless drawers promote simple cleaning for infection prevention
- Expansive work area with optional slide-out work surface
- 10-year limited warranty

Step 1

Choose a Base Cart Option:

Choose a cart, then add an accessory package to configure the cart for your specific requirements.

Procedure/Treatment Cart with Keylock



Standard
AM-PT-STD-KEY



Intermediate
AM-PT-INT-KEY



Compact
AM-PT-CMP-KEY

Procedure/Treatment Cart with Electronic Lock



Standard
AM-PT-STD-ELOK



Intermediate
AM-PT-INT-ELOK



Compact
AM-PT-CMP-ELOK

Cart Specifications

	Surface Height	Depth	Width
Standard	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
Intermediate	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
Compact	36" / 91 cm	24" / 61 cm	31" / 79 cm

Drawer Specifications

	Height	Depth	Width
Three Inch	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
Six Inch	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
Ten Inch	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

Specifications, availability & components are subject to change without notice.

Step 2

Choose an Accessory Package:

Choose your accessory package on the **next page**.

Step 2

Convenient Accessory Packages


Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

Avalo® Procedure/Treatment Cart Accessory Package 1

AM-PT-ACCPK1



3" Drawer Divider (1)



Writing Surface




Waste with Lid

Avalo® Procedure/Treatment Cart Accessory Package 2

AM-PT-ACCPK2



3" Drawer Divider (1)



Writing Surface



Waste with Lid



I.V. Pole



Sharps Container



Glove Box

Avalo® Procedure/Treatment Cart Accessory Package 3

AM-PT-ACCPK3



3" Drawer Divider (1)



Writing Surface



Articulating Arm



Laptop Security Tray

Can't find what you're looking for?

If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.