

Medication Carts

Avalo® Series Punch Card

Avalo Series medication carts excel in the areas that make the difference for nurses, like easy usage and access, a wide range of drawer configurations, enormous capacity, and smooth drawer action.

Available in standard light crème panels, or as Avalo Woodblend

Model PCS

43.75"(h) 24"(d) 31.5"(w)
Shown: Blush Salmon accent



Model PCL

43.75"(h) 24"(d) 38"(w)
Shown: Woodblend-Mountain Ash



Model PCXL

43.75"(h) 24"(d) 48.5"(w)
Shown: Woodblend-Northern Cherry



Offers a practical capacity of 350 punch cards plus room for liquids and supplies. Maximum capacity 450 punch cards.



Offers a practical capacity of 450 punch cards plus room for liquids and supplies. Maximum capacity 600 punch cards.



Offers a practical capacity of 600 punch cards plus room for liquids and supplies. Maximum capacity 750 punch cards.

Avalo Accent Colors



Taupe



Horizon Blue



Blush Salmon

Avalo Woodblend Finishes



Northern Cherry



Mountain Ash

Lock System Choices

Choose the lock system that best meets the environment's specific needs.

AutoLock™ Automatic Relocking System



Automatically secures a motionless unlocked medication cart with a user-defined relock time delay from 1-99 minutes. The system also permits the use of up to 5,000 individual user codes.

Keyless Access with Manual Relocking



Simply enter a 4-digit access code, push down the lock lever and the cart is ready for access. To secure the medication cart, push in the lock lever. Clear visual indication of lock status is apparent from a distance.

Key Access with Manual Relocking (standard)



Featuring the Best® Core Removable Lock System, this basic yet proven lock system provides your medication cart with the flexibility to change keys simply and easily. *Standard on Avalo carts.*

Ask about available security system upgrades

Accessories and Options

* Included on all Avalo Series Medication Carts. S & L Models also include Slide-Out Surface. XL Model includes Tracking Caster.

For additional accessories see CapsaHealthcare.com



Core Removable Lock*



Waste with Lid*




Narcotic Storage*



5" Brake Caster*



Cart Handle



Slide-Out Surface



Drawer Dividers



Sharps/Glove Box



Clear Top Mat



Storage Module



Punch Card Tub



Laptop Arm & Tray