

# AED WALL MOUNT CASES AND RACKS



## FEATURES

- Facilitates accessible and visible AED placement
- High tech design with large acrylic window provides a clear view of pad expiration date and the alarm activates when the AED is lifted out of the case
- RescueReady® indicator (assuring first time, every time shock delivery)



## AED VISIBILITY – EASY EMERGENCY ACCESS

### Introducing the NEW Cardiac Science Wall Mount Cases.

Because significantly reducing time to defibrillation saves lives from sudden cardiac arrest in a busy office or school environment, AED placement is critical. Cardiac Science has provided a solution with its innovatively designed wall mount cases and wire racks. Available in a variety of systems, the cases and racks facilitate accessible and visible AED placement.

The new Cardiac Science wall mount case is designed to hold the Cardiac Science AED and carry case, Ready Kit, and spare battery bag. This new high tech wall mount case is uniquely equipped with a strobe light alarm that is activated when the AED is lifted out of its cradle. The large, front paneled window allows users to easily check the defibrillation pad expiration date and RescueReady® indicator to ensure that the AED is ready for a rescue. The alarm strobe can be turned OFF with a key lock on the outside of the wall case, or by merely reinstalling the AED in the wall case cradle.

This wall mount is ideal for busy facilities such as corporations, airports, schools, health clubs, hotels, and manufacturing plants. Federal, state and municipal facilities are choosing Cardiac Science's wall mount cases and wire racks visibility, easy access, and innovative design. Creating Heartsafe Environments just got easier.

To learn more about Cardiac Science AED wall mount cases and wire racks, visit our website [www.cardiacscience.com](http://www.cardiacscience.com), or email [info@cardiacscience.com](mailto:info@cardiacscience.com), or call 1.866.289.5649 (USA and Canada only).









CARDIAC SCIENCE

CREATING HEARTSAFE ENVIRONMENTS

## ORDER INFORMATION

PART NUMBER DESCRIPTION

180-2021-001	<p><b>AED Wall Mount Storage Case with strobe light and alarm</b></p> <p>Exclusive, Cardiac Science AED wall mount storage case with strobe light and alarm. A superior storage wall case that perfectly holds the Cardiac Science AED with its carry case, Rescue Kit, and spare battery bag. The alarm strobe is only activated ON when the AED is lifted out of the wall case cradle. The alarm strobe can be turned OFF with a key lock on the outside of the wall case, or by merely reinstalling the AED in the wall case cradle. One key type per wall case is available, and two keys are included with the cabinet, chained inside the wall case. The wall case is surface mount. Wall case dimensions: 18.5" wide x 18" tall x 7.5" deep. Weight: 15.8 lbs.</p>	
180-2020-001	<p><b>AED Wall Mount Storage Case without strobe light and alarm</b></p> <p>AED wall mount storage case with no alarm or strobe light holds the AED for wall placement. This cutting edge, new design fits in a wide range of environments. The wall case is surface mount. Wall case dimensions: 18.5" wide x 18" tall x 7.5" deep. Weight: 15.8 lbs.</p>	
180-2021-002	<p><b>AED Wall Mount Storage Case with Security with strobe light and alarm</b></p> <p>AED wall mount storage case enhanced with security feature. When the AED is lifted out of the wall case cradle and the alarm strobe is activated ON, the pre-connected security features alerts emergency services. The wall case is surface mount. Wall case dimensions: 18.5" wide x 18" tall x 7.5" deep. Weight: 16.0 lbs.</p>	
180-2021-003	<p><b>Semi-Recessed AED Wall Mount Storage Case with strobe light and alarm</b></p> <p>AED Wall mount storage case embeds into wall. The semi-recessed feature makes the case suitable for a wide range of interior spaces and installation convenient. Wall case dimensions: 18.5" wide x 18" tall x 7.5" deep. Weight: 16.6 lbs.</p>	
170-2145-001	<p><b>Wall Mount Wire Rack for AED without straps</b></p> <p>AED wall mount wire rack tightly holds the AED for hanging on a wall. This wire rack provides convenient storage and easy accessibility of the AED during a rescue. Wire rack dimensions: 8.5" wide x 10.5" tall x 3.9" deep. Weight: 0.2 lbs.</p>	
170-2146-001	<p><b>Wall Mount Wire Rack for AED with straps</b></p> <p>AED wall mount wire rack tightly secures the AED in place with straps for accessible placement for hanging on a wall. Wall mount rack provides convenient storage and easy accessibility of the AED unit during a rescue. This wire rack with straps has been ambulance tested. Wire rack dimensions: 8.5" wide x 10.5" tall x 3.9" deep. Weight: 0.25 lbs.</p>	
170-2150-001	<p><b>Wall Mount Basket for AED without straps</b></p> <p>Basket securely holds the AED in place for hanging on a wall. Wall mount basket provides convenient storage and easy accessibility of the AED during a rescue. Wire rack dimensions: 14.25" wide x 9.3" tall x 6.8" deep. Weight: 0.8 lbs.</p>	
170-2152-001	<p><b>Wall Mount Basket for AED with straps</b></p> <p>Basket securely holds the AED in place for hanging on a wall. Wall mount basket provides convenient storage and easy accessibility of the AED during a rescue. This wall mount basket with straps has been ambulance tested. Wire rack dimensions: 14.25" wide x 9.3" tall x 6.8" deep. Weight: 0.85 lbs.</p>	

For more information contact Cardiac Science at:

**Corporate Headquarters**  
1900 Main Street, Suite 700, Irvine, CA, 92614, USA  
Phone: +1.949.797.3800 Fax: +1.949.951.7315

**Minnetonka Operations**  
5474 Felt! Road, Minnetonka, MN, 55343, USA  
Phone: +1.952.939.4181 Fax: +1.952.939.4191

**International Operations**  
Kirke Vaerloesevej 14, DK-3500 Vaerloese, Denmark  
Phone: +45.4438.0500 Fax: +45.4438.0501

**Cardiac Science Holdings (UK) Limited**  
Unit 13, Crown Industrial Estate, Canal Road, Timperley  
Altrincham, Cheshire, WA14 1TF, United Kingdom  
Phone: +44.161.926.0000 Fax: +44.161.926.0002

**Cardiac Science Deutschland GmbH**  
Max Planck Strasse 7B  
D-52249 Eschweiler, Germany  
Phone: +49.2403.5046.53 Fax: +49.2403.5046.54

**Worldwide Order Entry:**  
Toll Free: +1.800.991.5465 (USA and Canada)  
Phone: +1.949.797.3806  
Fax: +1.949.851.4418  
Email: aedorders@cardiacscience.com

**USA/Canada Customer Service/Technical Support:**  
Toll Free: +1.888.466.8686  
Phone: +1.949.797.3807  
Fax: +1.949.851.4425  
Email: aedsupport@cardiacscience.com

**International Customer Service/Technical Support:**  
Phone: +44.161.926.0000  
Fax: +44.161.926.0002  
Email: aed@cardiacscience.co.uk



CARDIAC SCIENCE

Nasdaq: DFIB

1.866.289.5649

WWW.CARDIACSCIENCE.COM

EMAIL: INFO@CARDIACSCIENCE.COM

