

CAIRE® Companion™ 5

5 LPM STATIONARY OXYGEN CONCENTRATOR

THE MOST RELIABLE AND TECHNOLOGICALLY ADVANCED

A high-performing stationary with an advanced design built specifically with longevity and durability in mind, the Companion 5™ is not only simple to use for patients, but its advanced technology offers a low total cost of ownership.

• Built To Last

Reduce costly breakdowns, after hours calls, and downtime with the unique design of the Companion 5™. It features a high-capacity compressor and high-capacity sieve beds equipped with proprietary moisture-resistant MGB SecoBlend™ sieve mixture. This unique design prolongs the life of the major components.

• autoFLOW™ Technology

Patent-pending autoFLOW™ technology (variable valve timing based on flow rate) was created specifically to help optimize performance and extend the life of the major components. It reduces the wear on the compressor by 25 percent at 2 LPM, decreases the operating temperatures to improve cooling, significantly reduces power draw, and improves the overall life of the concentrator.

• ECO Friendly

Drawing only 285 W at 2 LPM thanks to the reductions of autoFLOW™ Technology, the Companion 5 gives patients significant savings on their electrical costs, and helps to decrease wear on the electrical components.

• Ease of Service

With lengthy 2-year maintenance intervals, there is no gross particle filter to clean and no 9-volt battery to replace. The intake filter is easily accessible without opening the case through a filter door. The large and easy to read hour meter displays diagnostic alarm codes. These features help reduce after-hours calls and help over-the-phone diagnosis.




CAIRE®
A Chart Industries Company

2200 Airport Industrial Drive, Suite 500

Ball Ground, GA 30107

770.721.7759 | Toll Free 1.800.482.2473

www.CAIREMedical.com


CHART®

Copyright 2016 Chart Industries
ML-CONC0086 Rev A

PROVIDER FEATURES

- Advanced technology and design to minimize breakdowns
- autoFLOW™ Technology extends life of major components
- 2-year maintenance interval with filter door access
- Large hour meter displays diagnostic alarm codes
- Minimal internal components and modular design for easy service

PATIENT FEATURES

- Small profile fits easily into home
- Low power draw and electrical savings
- Easy-access humidifier stand
- Simple to read controls and alarms
- Low noise output with ultraSILENT technology
- No intake filter to clean

SAFETY FEATURES

- High flow alarm to prevent overdraw
- Low flow alarm to prevent kinked tubing
- Optional O₂ Monitor
- Flame-retardant plastic
- Compressor pressure relief valve
- Robust metal cannula barb and recessed flow



SPECIFICATIONS

Flow Rates	0.5 – 5.0 LPM
O₂ Concentration	90 % (+5.5 / -3 %)
Dimensions	21.5 x 12.5 x 13.5 inches (54.6 x 31.8 x 34.2 cm)
Weight	36 lbs (16.3 kg)
Sound Level	40 dB(A) @ 2.0 LPM
Operating Environment	41 to 104°F (4 to 40°C), Up to 95% humidity
Storage Environment	-4 to 158°F (-20 to 70°C), Up to 95% humidity
Power Consumption	285 W @ 2LPM, 350 W Maximum
Alarms	Power Failure, High Flow, Low Flow, General Malfunction, Low O ₂ Concentration (Optional)
Maintenance Schedule	Intake Filter (2-year replacement)
Maximum Tubing	50 ft (15.2 m)
Warranty	3 years

ORDERING INFORMATION

Companion 5™ Includes: Stationary Concentrator, Humidifier Adaptor Tubing, User Manual

15067005 120V / 60Hz, OM, USA Cord

15067013 120V / 60Hz, Non-OM, USA Cord

15067021 120V / 60Hz, OM, EU Cord

other countries' power cords may be available - contact customer service



2200 Airport Industrial Drive, Suite 500

Ball Ground, GA 30107

770.721.7759 | Toll Free 1.800.482.2473

www.CAIREMedical.com



Copyright 2016 Chart Industries
ML-CONC0086 Rev A