

Classic, Space Saver Table with Flat Top and Leg Rest



MODEL 1288-27

- No additional charge for flat top
- Classic wood-based construction
- Laminate hard surfaces with ultra-durable, PVC-overlay legs
- Flat, seamless rounded corner top
- Leg rest has 4 settings of approximately: 0°, 15°, 28° & 41°
- Full, laminate shelf
- 2" firm foam padding
- Heavyweight, knit-backed upholstery
- Easy assembly
- 350 lbs. load capacity under normal use (159 kg)

SPECIFICATIONS

1288-27	54"/72" Extended L x 27" W x 31" H
Top	54" L x 27" W
Leg Rest	18" L x 22" W

Prop 65 Warning



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Options



Laminate Colors



Dark Cherry



Maple



Gray

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black