

Taping Table with Shelf



MODEL 1702

- Heavy duty taping table with solid wood legs, H-braces and laminated full shelf in Natural
- 2" foam padding
- Natural finish
- 350 lbs. load capacity under normal use

SPECIFICATIONS

1702-27	48" L x 27" W x 36" H
1702-30	48" L x 30" W x 36" H

California residents please be advised of the following, pursuant to Proposition 65:

! **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Options



034- Paper Dispenser & Cutter Combo

040- Paper Cutter

030- Paper Dispenser

Finishes



Natural

Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black