


## Hi-Lo Mat Platform with Upholstered Top

### MODELS: 253-47/253-57/253-68

- Hi-Lo mat platform with easy-clean, ABS plastic, fully enclosed base
- Electric power height adjustment with foot control
- Ultra-strong, uniframe base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Custom bumper strip helps protect top's edges from tears
- 2" ultra-firm padding
- Low height may help achieve ADA tax credit
- Made in the USA by Clinton Industries
- 600 lbs. load capacity under normal use



#### MODEL 253-47 SPECIFICATIONS

Top	7' L x 4' W x 22"—34" H
-----	-------------------------

#### MODEL 253-57 SPECIFICATIONS

Top	7' L x 5' W x 22"—34" H
-----	-------------------------

#### MODEL 253-68 SPECIFICATIONS

Top	8' L x 6' W x 22"—34" H
-----	-------------------------

**IMPORTANT NOTE:** Due to the size of the product, it **MUST** ship to a truck height dock.



California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

### Standard Upholstery



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black