

Hi-Lo, Hydraulic, Mat Platform with Upholstered Top

MODELS: 270-46/270-47/270-57/270-68

- Ultra-strong, powder-coated steel frame
- Hydraulic height adjustment
- Hands-free, electric, foot control
- Custom bumper strip helps protect top's edges from tears
- Models 270-46 & 270-47 have seamless upholstery
- Models 270-57 & 270-68, because of their size, have upholstery that is seamed in the middle
- 2" ultra-firm padding
- Low height may help achieve ADA tax credit
- Made in the USA by Clinton Industries
- 600 lbs. load capacity under normal use

IMPORTANT NOTE: Due to the size of the 270-68, it **MUST** ship to a truck height dock.



MODEL 270-46 SPECIFICATIONS

Top 6' L x 4' W x 20"—30" H

MODEL 270-47 SPECIFICATIONS

Top 7' L x 4' W x 20"—30" H



MODEL 270-57 SPECIFICATIONS


Top 7' L x 5' W x 20"—30" H

MODEL 270-68 SPECIFICATIONS

Top 8' L x 6' W x 20"—30" H



California residents please be advised of the following, pursuant to Proposition 65:

1
 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Upholstery Colors



Black



Burgundy



Country Mist



Desert Tan



Gunmetal



Royal Blue



Slate Blue



Warm Gray



Wedgewood