

Platform Mounted Parallel Bars

MODELS: 3-2007 / 3-2010 / 3-2012


- Parallel bars mounted to 1-1/2" thick hardwood plywood base with ultra-durable, UV-cured satin finish
- Safety tread at both ends of the platform
- Width adjustable, 1-1/2" one-piece, stainless steel handrails
- Heavy duty, square, telescoping, steel uprights
- Easy height changes with spring-loaded plungers and numbered height indicator strips
- Second height knob removes side to side play
- End bumpers on all handrails
- 400 lbs. load capacity under normal physical therapy use



IMPORTANT NOTE: *Intended for physical therapy use only, not for gymnastics use.*



California residents please be advised of the following, pursuant to Proposition 65:

1
 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

MODEL 3-2007 SPECIFICATIONS

7' L x 18"-28" W range x 26"-39" H range

MODEL 3-2010 SPECIFICATIONS

10' L x 18"-28" W range x 26"-39" H range

MODEL 3-2012 SPECIFICATIONS

12' L x 18"-28" W range x 26"-39" H range