

Heavy-Duty, Parallel Bars

MODEL 3-2106

- Parallel bars are mounted to a 1-1/2" thick hardwood plywood base with ultra-durable UV-coated satin finish
- Safety tread at both ends of the platform
- Middle upright with through pin lock for extra support
- Extra-wide 48" width for bariatric patients
- Width adjustable 1-1/2" one-piece, stainless steel handrails
- Easy height adjustment with spring loaded plungers and numbered height strip
- Heavy duty, square telescoping steel uprights
- Second height knob removes side to side play
- End bumpers on all handrails
- 500 lbs. load capacity under normal physical therapy use



IMPORTANT NOTE: *Intended for physical therapy use only, not for gymnastics use.*



California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

SPECIFICATIONS

10' L x 20"-34" W range x 26"-39" H range