

Wall Mounted Folding Parallel Bars

MODEL 3-3307

- Stainless steel handrails
- Requires only 14" of floor space when stored
- Folds up against the wall when not in use (folded height 94"). Please allow for a minimum 8' ceiling height
- 1-1/2" one-piece, stainless steel handrails
- Easy height adjustment with spring loaded plungers and numbered height stripon front legs
- Heavy duty, square telescoping steel front legs and stainless steel wall uprights
- Second height knob removes side to side play in front legs
- End bumpers on all handrails
- 400 lbs. load capacity under normal physical therapy use

IMPORTANT NOTE: *Intended for physical therapy use only, not for gymnastics use.*

NOTE: *Mounting hardware not included. Check with your building manager for the proper type of fasteners needed for your specific wall structure.*

SPECIFICATIONS

Deployed	7' L x 28" W x 28"-42" H range
Folded	94" High



California residents please be advised of the following, pursuant to Proposition 65:

! **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



Deployed



Stored