

Folding Parallel Bars

MODELS: 3-3317 / 3-3310

- Width adjustable floor folding parallel bars with double stiffening rails
- Width adjustable, 1-1/2" one-piece, stainless steel handrails
- Easy height adjustment with spring loaded plungers and numbered height strip
- Heavy duty, square telescoping steel uprights
- Second height knob removes side to side play
- End bumpers on all handrails
- 350 lbs. load capacity under normal physical therapy use

IMPORTANT NOTE: *Intended for physical therapy use only, not for gymnastics use.*



California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

MODEL 3-3317 SPECIFICATIONS

7' L x 16"-24" W range x 28"-41" H range

MODEL 3-3310 SPECIFICATIONS

10' L x 16"-24" W range x 28"-41" H range