


PARALLEL BARS

**Wall Mounted Folding Parallel Bars**

- Stainless steel handrails
- Requires only 14' of floor space when in use
- Folds up against the wall when not in use
- Folded Height 7" 10" (238.76 cm)
- 400 lbs. load capacity (181.43 kg)



3-3307	length	width	height range
	7'	28"	28" – 42"
	2.13 m	71.12 cm	71.12 cm – 106.68

**Floor Mounted Parallel Bars**

- Stainless steel adjustable handrails
- Same features as platform-mounted model, except it is floor-mounted
- 400 lbs. load capacity when properly mounted (181.43 kg)

3-4007	length	width range	height range
	7'	adjusts 6"/side	26"–39"
	2.13 m	adjusts 15.23 cm/side	66 cm – 99 cm

3-3307



3-4007

3-4010	10' long model (3 m)
--------	----------------------

3-4012	12' long model (3.65 m)
--------	-------------------------

**Bariatric Parallel Bars**

- Height and width adjustment
- 3rd upright with through pin lock for extra support
- Extra width for bariatric patients
- 48" platform width (121.92 cm)
- 500 lbs. load capacity (226.79 kg)

3-2106	length	width range	height range
	10'	20" – 34"	26"–39"
	3.04 m	50.80 cm – 86.36 cm	66 cm – 99 cm



3-2106

**Folding Parallel Bars**

- Stainless steel handrails
- Unique double stabilizer bars stiffen the rails
- Spring-loaded plungers and sure-lock positioning knobs
- 350 lbs. load capacity under normal use (158.57 kg)

3-3317	length	width range	height range
	7'	16" – 24"	28"–41"
	2 m	41 cm – 61 cm	71 cm – 104cm

3-3310	10' long model (3.04 m)
--------	-------------------------



3-3317



Folded

**Optional Center Obduction Board**

- Available with platform style only
- Easily removes when not needed

045



045