

Floor Mounted Parallel Bars

MODELS: 3-4007 / 3-4010 / 3-4012

- Parallel bars and uprights designed to mount directly to the floor
- Width adjustable, 1-1/2" one-piece, stainless steel handrails
- Heavy duty, square, telescoping, steel uprights
- Easy height changes with spring-loaded plungers and numbered height indicator strips
- Second height knob removes side to side play
- End bumpers on all handrails
- 400 lbs. load capacity under normal physical therapy use when properly installed

IMPORTANT NOTE: *Intended for physical therapy use only, not for gymnastics use.*



California residents please be advised of the following, pursuant to Proposition 65:

! **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

MODEL 3-4007 SPECIFICATIONS

7' L x adjusts 6"/side W x 26"-39" H range

MODEL 3-4010 SPECIFICATIONS

10' L x adjusts 6"/side W x 26"-39" H range

MODEL 3-4012 SPECIFICATIONS

12' L x adjusts 6"/side W x 26"-39" H range