

## BOLSTERS/WEDGES/PILLOWS

### Bolster

50	diameter	length
	6" 15.24 cm	24" 60.96 cm

### Bolster

51	diameter	length
	8" 20.32 cm	24" 60.96 cm

### Half Round with Flat Side

52	length	width	height
	18" 46.0 cm	8" 20.32 cm	6" 15.24 cm

### Cube

53	length	width	height
	12" 30.5 cm	8" 20.32 cm	3" 7.62 cm

### Cube

54	length	width	height
	12" 30.5 cm	8" 20.32 cm	6" 15.24 cm

### Wedge

56	length	width	height
	18" 46.0 cm	24" 60.96 cm	6" 15.24 cm

### Wedge

57	length	width	height
	32" 81.28 cm	20" 50.8 cm	12" 30.48 cm

### Wedge

58	length	width	height
	16" 40.64 cm	8" 20.32 cm	8" 20.32 cm

### Cube/Incline

59	length	width	height
	32" 81.28 cm	10" 25.4 cm	12" 30.48 cm

### Cube/Incline

61	length	width	height
	32" 81.28 cm	20" 50.8 cm	12" 30.48 cm

### Large Pillow


60	length	width	height
	13" 33 cm	20" 50.8 cm	3" 7.62 cm

### Small Pillow

20	length	width	height
	14" 35.56 cm	12" 30.48 cm	3" 7.62 cm

Clinton bolsters, wedges and pillows feature:

- Firm high density urethane foam core
- Knit-backed, sewn, vinyl cover in choice of color
- Easy-clean surface



## VANGUARD SERIES/RACS

### Astoria MaxRac with Mirror

- 6 independent, quick-change positions for holding your choice of exercise band (50 yds, (45.72 m) hollow core center only). All items shown not included.
- Minimum floor space requirement
- Large interior space with door
- Mirror has polished edges & ANSI 297 safety backing
- Raised edge on top provides extra storage space
- Cuff weight "sliding hook system" makes lateral positioning easy
- Strong polycarbonate dumbbell rods are angled for added stability
- Protective bumper moulding around base
- 4" swivel casters/2 locking (10.16 cm)
- Holds up to 64 cuff weights, 22 dumbbells and 6 rolls of exercise band

5101B	width	depth	height
	28" 71.28 cm	28" 71.28 cm	65" 165.1 cm

5101A Same as 5101B without Mirror



Back



Raised Edge on Top

### Bayside MaxRac with Mirror

- 6 independent, quick-change positions for holding your choice of exercise band (50 yds, (45.72 m), hollow core center only). All items shown not included.
- Minimum floor space requirement
- Large interior space with door
- Mirror has polished edges & ANSI 297 safety backing
- Raised edge on top allows for extra storage space
- Cuff weight "sliding hook system" makes lateral positioning easy
- Strong polycarbonate dumbbell rods are angled for added stability
- Protective bumper moulding around base
- 4" swivel casters/2 locking (10.16 cm)
- Holds up to 48 cuff weights, 22 dumbbells 20 weight bars and 6 rolls of exercise band

5102B	width	depth	height
	28" 71.28 cm	28" 71.28 cm	65" 165.1 cm

5102A Same as 5102B without Mirror



Back



LAMINATE COLOR CHOICES

Gray—1GR

Natural—1NT

Maple—1MP

Dark Cherry—1DC

White—1WH