

## SC Series, Extra-Wide, Padded, Blood Drawing Chair with Padded Flip Arms

### MODEL 6006-U

- Contoured, padded seat and 2" padded backrest
- Height and depth adjustable, dual, padded, flip arms
- Flip arm positions in place when in use, then flips out of the way for easy exit
- Dual ClintonClean™ elbow rests
- Contemporary styling
- Gray powder-coated frame finish
- Adjustable height and depth flip arm
- Non-marking, nylon glides
- 300 lbs. load capacity under normal use (136 kg)

**IMPORTANT NOTE:** Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

**ORDERING NOTE:** MUST add arm suffix number to model number when ordering with 665 arm option, otherwise the chair will be shipped as shown.

**All blood drawing products are non-returnable unless they are unused in unopened, original packaging.**



#### SPECIFICATIONS

Seat	30" W x 18" D x 18" H
Overall	45" W x 25" D x 32.5"H
Arm Height Range	27" — 33"
Arm Adjustment	6"



California residents please be advised of the following, pursuant to Proposition 65:

**1**  
**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

#### OPTIONS



665 ClintonClean™, Armrests and Flip Arm(s)

#### Standard Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black