

CabinetRac with Doors and Drop Leaf

- 6 independent, quick-change positions for holding your choice of exercise band (50 yds. (45.72 m), hollow core center only). All items shown not included.)
- Drop leaf feature is ideal for writing/support surface
- 2 locking doors with 3 adjustable shelves
- Storage nook for Velcro® with door
- Mirror has polished edges and ANSI 297 safety backing
- Strong polycarbonate dumbbell rods are angled for added stability
- Protective bumper moulding around edge
- 4" swivel casters/2 locking (10.16 cm)
- Holds up to 48 cuff weights, 22 dumbbells, 4 rolls of Velcro®, 6 rolls of exercise band and much more

6150A	width	depth	height
	30"	30"	73"
	76.2 cm	76.2 cm	185.42 cm



MirrorRac

- 6 independent, quick-change positions for holding your choice of exercise band (50 yds. (45.72 m), hollow core center only). All items shown not included.
- Mirror with ANSI 297 safety backing
- Strong polycarbonate dumbbell rods are angled for added stability
- Protective bumper moulding around edge
- 4" swivel casters/2 locking (10.16 cm)
- White pegboard to hold cuff weights
- Holds up to 48 cuff weights, 22 dumbbells, 6 rolls of exercise band and much more (Not included)

6154	width	depth	height
	31 1/2"	18"	72"
	80 cm	45.72 cm	182.88 cm



CabinetRac with Doors Plus

- Versatile, multi-item weight/storage cabinet with 2 locking doors
- 6 independent, quick-change positions for holding your choice of exercise band (50 yds. (45.72 m), hollow core center only). All items shown not included.)
- 6 shelves (3 adjustable and 3 fixed)
- Open top storage nook
- Strong polycarbonate dumbbell rods are angled for added stability
- Protective bumper moulding around edge
- 4" swivel casters/2 locking (10.16 cm)
- White pegboard to hold cuff weights
- Holds up to 48 cuff weights, 22 dumbbells, 6 rolls of exercise band and much more

6151	width	depth	height
	29"	23"	72"
	73.66 cm	58.42 cm	182.88 cm



CabinetRac with Doors

- Versatile, multi-item weight/storage cabinet with 2 locking doors
- 6 independent, quick-change positions for holding your choice of exercise band (50 yds. (45.72 m), hollow core center only). All items shown not included.
- 4 shelves (3 adjustable, 1 fixed)
- Strong polycarbonate dumbbell rods are angled for added stability
- Protective bumper moulding around edge
- 4" swivel casters/2 locking (10.16 cm)
- White pegboard to hold cuff weights
- Holds up to 48 cuff weights, 22 dumbbells, 6 rolls of exercise band and much more (Not included)

6155	width	depth	height
	29"	23"	62"
	73.66 cm	58.42 cm	157.48 cm



LAMINATE COLOR CHOICES



LAMINATE COLOR CHOICES

