

## Series 600, Low Height, Classic, Treatment Table



### MODEL 6196

- Plywood top and solid hardwood legs
- 6-leg design with triple bolted corner legs
- Built-in middle leg levelers
- Steel reinforced frame
- Laminate full length shelf with H-braces
- Pillow included
- 2" firm foam padding
- Lower height makes patient access easier
- 600 lbs. capacity under normal use



### SPECIFICATIONS

6196

78" L x 30" W x 26" H

California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

### Options



034- Paper  
Dispenser & Cutter  
Combo

040- Paper Cutter

030- Paper  
Dispenser

### Finishes



Natural



Dark Cherry

### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black