

Mobile Adult Mirror with Cuff Weight Rac

MODEL 6210-9

- 1/4" thick distortion-free glass mirror with ANSI 297 safety backing
- All laminate construction
- Rolls easily on 3" swivel casters
- ABS plastic bumpers on base sides
- Cuff weight sliding hook system makes lateral positioning fast and easy. 12 cuff weight hooks included
- Holds up to 48 cuff weights (Not included)



California residents please be advised of the following, pursuant to Proposition 65:

Prop 65 Warning

! **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



MODEL 6210-9 SPECIFICATIONS

Overall	31.5" W x 20.5" D x 75.75" H
Base	31.5" W x 20.5" D
Top	28.25" W x 4" D x 70.75" H

Laminate



Maple