

## Lab X Series, Blood Drawing Chair with Padded Flip Arm and Drawer

### MODEL 66029-P

- E-Z-Clean, one-piece, black, plastic seat
- Height and depth adjustable, ClintonClean™ flip arm
- Flip arm positions in place when in use, then flips out of the way for easy exit
- ClintonClean™ elbow rest
- ClintonClean™ drawer top
- Height-adjustable, one-piece, plastic drawer
- 1-1/4" square, heavy duty, all-welded, tubular steel, frame (3.18 cm)
- Neutral Gray powder-coated frame
- 300 lbs. load capacity under normal use (136 kg)

**NOTE:** If padded upholstered arms are selected, use color chart to choose color.

**IMPORTANT NOTE:** Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

**ORDERING NOTE:** MUST add arm suffix number to model number when ordering with 600 arm option, otherwise the chair will be shipped as shown.

**All blood drawing products are non-returnable unless they are unused in unopened, original packaging.**



#### MODEL 66029-P SPECIFICATIONS

Seat	18" W x 17" D x 20" H
Overall	40" W x 27" D x 36"H
Arm Height Range	27" — 33"
Arm Adjustment	6"



California residents please be advised of the following, pursuant to Proposition 65:

**1**  
**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

#### OPTIONS



600 Padded, Armrest and Flip Arm

#### Standard Upholstery Colors for Optional Arms

