

## Lab X Series, Extra-Wide /Extra-Tall, Blood Chair/ Padded Flip Arm and Drawer

### MODEL 66092

- Padded seat and backrest
- Height and depth adjustable, padded, flip arm
- Flip arm(s) positions in place when in use, then flips out of the way for easy exit
- ClintonClean™ elbow rest
- ClintonClean™ drawer top
- Height-adjustable, one-piece, all-plastic drawer
- Extra height reduces bending for attendant
- A full 35" wide seat to accommodate a larger patient, or parent and child (89 cm)
- 1-1/4" square, heavy duty, all-welded, tubular steel frame (3.18 cm)
- Comfortable 2" padding throughout (5 cm)
- Neutral Gray, powder-coated frame
- Easy-clean, durable vinyl in a selection of colors
- ClintonClean™ arms option available
- 400 lbs. combined load capacity under normal use (181.4 kg)



**IMPORTANT NOTE: Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.**

**ORDERING NOTE: MUST add arm suffix number to model number when ordering any model with arm option, otherwise chair will be shipped as shown.**

**All blood drawing products are non-returnable unless they are unused in unopened, original packaging.**



California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

### MODEL 66092 SPECIFICATIONS

Seat	35" W x 18" D x 26" H
Overall	53" W x 26" D x 42" H
Arm Height Range	33" — 39"
Arm Adjustment	6"

### OPTIONS



665 ClintonClean™, Armrest and Flip Arm(s)

### Standard Upholstery Colors

