



# WHY CLINTON?

Buyers of medical equipment are faced with evaluating clinical needs and justifying expenditures on nearly every purchase. Now, Clinton makes it a whole lot easier on everybody with products that maximize performance and aesthetics while minimizing costs.

**Performance.** With nearly forty years in medical manufacturing, Clinton has met the evolving challenges of the industry with products that fit the needs of providers, while at the same time maintaining a track record of durability that boasts the best warranty in the business.

**Aesthetics.** Nobody likes outdated anything, so Clinton has changed its products over the years to reflect the advances in medical practice, design and contemporary decor, making products with appeal in both form and function.

**Cost.** Bells and whistles are great, but they also drive up costs. On the other hand, why pay more if you don't have to? Clinton products fill needs instead of empty pockets, providing medical professionals with just the right products—at just the right prices.

**IMPORTANT NOTE TO DISTRIBUTORS DOING BUSINESS IN CALIFORNIA:**

Clinton Industries products drop shipping or to our knowledge being reshipped into California will carry appropriate Safe Harbor Warnings in compliance with California Proposition 65.

For more information on California Proposition 65, go to: [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)



CLINTON  
**STORE & GO™  
CARTS**



Add more function and value to everyday task carts with a Clinton Store & Go Cart. Clinton's line of Store & Go Carts are not only practical for clinical tasks like phlebotomy, but also feature a wide range of storage options. So, one cart can cover a variety of daily assignments without stopping to resupply. Select from six, affordable models to suit your individual needs for trays, bins, compartments, drawers and security—all backed up with Clinton's 5-year limited warranty.

## 67001 Store & Go Cart

- ClintonClean™ solid plastic top with adjustable tray brackets
- Easy-grip, ClintonClean™ solid plastic handle
- One, storage tray with two, plastic, storage bins
- Stable, 5-leg, chrome base
- Knob, height adjustable, chrome post
- White, powder-coated steel, glove box holder
- Universal attachment rail
- 3" dual-wheel, hooded casters
- Fast, easy assembly with basic tools



Integrated handle and adjustable tray brackets



Tray pictured on top and all supplies shown not included, and not sold by Clinton.

### 67001

Top: 20.5" x 16"  
Top brackets can accommodate trays from 16" x 15" up to 17.75" x 12"  
Small Bin Holders: 7.5" x 5" on each side of the post  
Small Plastic Bins (2): 6.75" x 3.5"  
Base: 5-leg, 22.5" diameter  
Height Range: 26.5" x 37.25"

## 67021 Store & Go Cart

- ClintonClean™ solid plastic top with adjustable tray brackets
- Easy-grip, ClintonClean™ solid plastic handle
- Two, mid-height storage compartments featuring: one, storage tray with two, plastic, storage bins and a drop-in, large, plastic bin
- Stable, 5-leg, chrome base
- Knob, height adjustable, chrome post
- White, powder-coated steel, glove box holder
- Universal attachment rail
- 3" dual-wheel, hooded casters
- Fast, easy assembly with basic tools



Integrated handle and adjustable tray brackets



Tray pictured on top and all supplies shown not included, and not sold by Clinton.

### 67021

Top: 20.5" x 16"  
Top brackets can accommodate trays from 16" x 15" up to 17.75" x 12"  
Small Bin Holders: 7.5" x 5" on each side of the post  
Small Plastic Bins (2): 6.75" x 3.5"  
Medium Bin holder 10" x 7.25"  
Medium plastic bin: 9.75" x 7.125"  
Base: 5-leg, 22.5" diameter  
Height Range: 26.5" x 37.25"

All blood drawing products are non-returnable, unless they are unused in unopened, original packaging.

## 67002 Store & Go Cart

- Easy-clean, moulded plastic, rimmed top with convenient plastic storage drawer
- Easy-grip, ClintonClean™ solid plastic handle
- Two, mid-height, storage compartments featuring: one, storage tray with two, plastic, storage bins and a drop-in, large, plastic bin
- Large, stable, 5-leg, chrome base
- Pneumatic, height adjustable cylinder
- White, powder-coated, steel, glove box holder
- Universal attachment rail
- 4" dual-wheel, hooded casters
- Fast, easy assembly with basic tools



ClintonClean™ handle, and Easy-clean, plastic drawer

## 67031 Store & Go Cart

- ClintonClean™ solid plastic top with adjustable tray brackets
- Easy-grip, ClintonClean™ solid plastic handle
- Two, mid-height storage compartments featuring: one, storage tray with two, plastic, storage bins, a large, tray section with a large, plastic bin and a lockable\* storage bin
- Large, stable, 5-leg, chrome base
- Pneumatic, height adjustable cylinder
- White, powder-coated, steel, glove box holder
- Universal attachment rail
- 4" dual-wheel, hooded casters
- Fast, easy assembly with basic tools



Integrated handle and adjustable tray brackets

Locking storage bin

All blood drawing products are non-returnable, unless they are unused in unopened, original packaging.



Tray pictured on top and all supplies shown not included, and not sold by Clinton.

### 67002

Top: 20.75" x 20" / Including Handle: 24" x 20"  
Small Bin Holders: 7.5" x 5" on each side of the post  
Small Plastic Bins (2): 6.75" x 3.5"  
Medium Bin Holder 10" x 7.25"  
Medium Plastic Bin: 9.75" x 7.125"  
Plastic Drawer: 17.25" x 7"  
Base: 5-leg, 26.5" diameter  
Height Range: 31.75" x 42"



Tray pictured on top and all supplies shown not included, and not sold by Clinton.

### 67031

Top: 20.5" x 16"  
Top brackets can accommodate trays from 16" x 15" up to 17.75" x 12"  
Small Bin Holders: 8.75" x 4" on each side of the post  
Small Plastic Bins (2): 6.75" x 3.5"  
Large Bin Tray 18.75" x 10.75"  
Large Plastic Bin: 15.5" x 9.75"  
Lockable Bin: 15.5" x 9.75"  
Base: 5-leg, 26.5" diameter  
Height Range: 31.125" x 41.375"

## 67022 Store & Go Cart

- Easy-clean, moulded plastic, rimmed top with convenient plastic storage drawer
- Easy-grip, ClintonClean™ solid plastic handle
- Two, mid-height storage compartments featuring: one, storage tray with two, plastic, storage bins, a large, tray section with a large, plastic bin and a lockable\* storage bin
- Large, stable, 5-leg, chrome base
- Pneumatic, height adjustable cylinder
- White, powder-coated steel, glove box holder
- Universal attachment rail
- 4" dual-wheel, hooded casters
- Fast, easy assembly with basic tools



ClintonClean™ handle



Plastic drawer



Lockable storage bin



Tray pictured on top and all supplies shown not included, and not sold by Clinton.

### 67022

Top: 20.75" x 20" / Including Handle: 24" x 20"  
 Small Bin Holders: 8.75" x 4" on each side of the post  
 Small Plastic Bins (2): 6.75" x 3.5"  
 Large Bin Tray 18.75" x 10.75"  
 Large Plastic Bin: 15.5" x 9.75"  
 Lockable Bin: 15.5" x 9.75"  
 Plastic Drawer: 17.25" x 7"  
 Base: 5-leg, 26.5" diameter  
 Height Range: 31.75" x 42"

All blood drawing products are non-returnable, unless they are unused in unopened, original packaging.

## 67032 Slimline Store & Go Cart

- ClintonClean™ solid plastic top with adjustable tray brackets
- Easy-grip, ClintonClean™ solid plastic handle
- Three, storage shelves
- Three, convenient, plastic drawers
- Large, locking storage compartment
- Stable, 4-caster, ClintonClean™ solid plastic base
- Top raises with a spring-loaded plunger, height adjustment
- Double, white, powder-coated steel, glove box holder
- 4" dual-wheel, hooded casters
- Fast, easy assembly with basic tools



Integrated handle with adjustable tray brackets



Three, easy-clean, plastic drawers



Lockable storage compartment with plastic bin



Tray pictured on top and all supplies shown not included, and not sold by Clinton.



Spring-loaded, plunger, height adjustment with 8 settings up to 43"

### 67032

Top: 20.5" x 16"  
 Top brackets can accommodate trays from 16" x 15" up to 17.75" x 12"  
 Top of Cart: 18" x 8.5"  
 Storage Shelves (3): 11" x 3.75"  
 Plastic Drawer (3): 17.25" x 7"  
 Lockable Storage Compartment: 15" x 8.25"  
 Plastic Storage Bin Inside Lockable Storage Compartment: 10.5" x 6"  
 Base: 22.5" x 14"  
 Height Range: (5 settings) 36", 37.5", 39", 40.5" & 42"

All blood drawing products are non-returnable, unless they are unused in unopened, original packaging.



525 East Market St., York, PA 17403

800-441-9131 • 717-848-3519

Toll-Free Fax 866-522-1233

Fax 717-843-5871

[www.clinton-ind.com](http://www.clinton-ind.com)