

67001 Store & Go Phlebotomy Cart

MODEL 67001

- ClintonClean™ solid plastic top with adjustable tray brackets
- Easy-grip, ClintonClean™ solid plastic handle
- One, storage tray with two, plastic, storage bins
- Stable, 5-leg, chrome base
- Knob, height adjustable, chrome post
- White, powder-coated steel, glove box holder
- Universal attachment rail
- 3" dual-wheel, hooded casters
- Fast, easy assembly with basic tools

NOTE: *Items pictured on top not included*

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.

SPECIFICATIONS

Top: 20.5" x 16"

Top brackets can accommodate trays from 16" x 15" up to 17.75" x 12"

Small Bin Holders: 7.5" x 5" on each side of the post

Small Plastic Bins (2): 6.75" x 3.5"

Base: 5-leg, 22.5" diameter

Height Range: 26.5" x 37.25"



California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1