

Easy to handle and maneuver in tight spaces with ample space for supplies and equipment, Clinton's phlebotomy carts are the right fit for busy nurses and phlebotomists.

5-Leg, Phlebotomy Carts feature:

- Gray laminate cabinet construction
- Easy-clean, vacuum-formed, gray, ABS, plastic top with textured finish and raised rim on all sides
- Lockable door with soft-close hinges
- Removable plastic storage tray(s) on easy slide tracks
- Easy grip handle 1" diameter, powder-coated, push handle
- Pneumatic height adjustment
- Stable, 5-leg, aluminum base
- Smooth rolling, 3" ball bearing casters (7.62 cm)
- Optional glove box holder available



One-Bin Phlebotomy Cart

67100			
top	width 20" 50.8 cm	depth 20 ³ / ₄ " 52.7 cm	
cabinet	width 13 ¹ / ₄ " 33.65 cm	depth 17 ¹ / ₄ " 43.81 cm	height 7 ³ / ₄ " 19.68 cm
base diameter	24" 60.96 cm		
height range	33" – 43" 83.82 cm – 109.22 cm		

Two-Bin Phlebotomy Cart

67200			
top	width 20" 50.8 cm	depth 20 ³ / ₄ " 52.7 cm	
cabinet	width 13 ¹ / ₄ " 33.65 cm	depth 18" 45.72 cm	height 10 ¹ / ₂ " 26.67 cm
base diameter	24" 60.96 cm		
height range	33" – 43" 83.82 cm – 109.22 cm		

H-Base, Phlebotomy Carts feature:

- Gray laminate cabinet construction
- Easy-clean, vacuum-formed, gray, ABS, plastic top with textured finish and raised rim on all sides
- Lockable door with soft-close hinges
- Removable plastic storage tray(s) on easy-slide tracks
- Easy grip, 1" powder-coated steel tubing handle
- Pneumatic height adjustment
- H-leg, powder-coated steel base
- Smooth rolling, 3" ballbearing casters (7.62 cm)
- Optional glove box holder available

Small, H-Base, Pneumatic, Two-Bin Phlebotomy Cart

67202			
top	width 20" 50.8 cm	depth 20 ¹ / ₂ " 52 cm	
cabinet	width 13 ¹ / ₄ " 33.65 cm	depth 18" 45.72 cm	height 10 ¹ / ₂ " 26.67 cm
height range	29" – 42 ³ / ₄ " 73.66 cm – 108.58 cm		

Large, H-Base, Pneumatic, Four-Bin Phlebotomy Cart

67203			
top	width 20" 50.8 cm	depth 20 ¹ / ₂ " 52 cm	
cabinet	width 13 ¹ / ₄ " 33.65 cm	depth 17 ³ / ₄ " 45 cm	height 19 ¹ / ₂ " 49.53 cm
height range	29" – 42 ³ / ₄ " 73.66 cm – 108.58 cm		



Two plastic bins slide on plastic channels



1" powder-coated, steel tubing handle

(Tray on top and supplies shown not included)



Four plastic bins slide on plastic channels

(Tray on top and supplies shown not included)



Optional glove box holder mounts to back (Order Single White Steel Glove Box Holder GW-2000)



One or two plastic bins slide on plastic channels

(Tray on top and supplies shown not included)