

Two-Bin Phlebotomy Cart

MODEL 67200

- Gray laminate cabinet construction
- Easy-clean, vacuum-formed, gray, ABS, plastic top with textured finish and raised rim on all sides
- Lockable door with soft close hinges
- 2 Removable plastic storage bins on easy slide tracks 16-1/2" L x 10-7/8" W x 3-1/2" D (42cm L x 28cm W x 9cm D)
- Easy grip, ClintonClean™ solid plastic handle
- Pneumatic height adjustment
- Stable, 5-leg, black nylon base
- Smooth rolling, 3" ball bearing casters (7.62 cm)
- Optional glove box holder available. Add product GW-2000 Single White Steel Glove Box Holder



NOTE: Items pictured on top not included

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.

Lock is not engineered to secure narcotics or other limited access drugs.

MODEL 67200 SPECIFICATIONS

Top	20" W x 20" D
Base Diameter	24"
Height Range	33" — 43"
Cabinet	13.5" W x 18" D x 10.5"


OPTIONS



GW-2000 Single White Steel Glove Box Holder



Prop 65 Warning

1
 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Laminate Color

Gray