

## Cool Pals, 2-Bin, Phlebotomy Cart

### MODEL 67231

- Clinton Kids™ Cool Pals graphics embedded in laminate will not tear or peel
- Lockable door on soft close hinges
- Two, 3.5" x 11" x 16.5" removable plastic storage bins on easy slide tracks (9 cm x 28 cm x 42 cm)
- Easy-clean, vacuum-formed, gray, ABS, plastic top with textured finish and raised rim on all sides
- Easy grip, ClintonClean™ solid plastic handle
- Pneumatic height adjustment
- Black, nylon 5-leg base
- Smooth rolling, 3" ball bearing casters (7.62 cm)
- Optional glove box holder available. Add product GW-2000 Single White Steel Glove Box Holder



**NOTE: Items pictured on top not included**

**All blood drawing products are non-returnable unless they are unused in unopened, original packaging.**

#### MODEL 67231 SPECIFICATIONS

Top	20" W x 20" D
Base Diameter	24"
Height Range	33" — 43"
Cabinet	13.5" W x 18" D x 10.5"



#### Prop 65 Warning

**!** **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

#### OPTIONS



GW-2000 Single White Steel Glove Box Holder