

Recliner Series Hi-Lo Blood Drawing Chair

MODEL 6810

- Smooth, fully electric reclining and height adjustment
- Self leveling arms and synchronized footrest that move with chair's seat and backrest
- Keeps patients comfortable and may better accommodate patients who may faint during a procedure
- Hand controller
- All-welded, steel, neutral gray powder-coated frame
- Heavy duty, all-welded, steel frame
- Padded seat, backrest, armrests and footrest
- Armrests flip up out of the way for easy entry and exit
- Adjustable foot levelers
- Easy-clean, durable, stain-resistant vinyl
- 375 lbs. load capacity under normal use (170 kg)

IMPORTANT NOTE: Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.



MODEL 6810 SPECIFICATIONS

Seat	20.5" W x 17" D x 22" – 30" H range
Back	20" W x 32" H with 25° – 75° angle
Footrest	19.5" W x 14" D with 5° – 90° angle
Overall	35" W x 50" D x 35" retracted/66" extended
Arm Adjustment	21" between arms



California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Standard Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black