

PHLEBOTOMY STANDS & WORK STATIONS

Phlebotomy stands and work stations are a quick, easy and economical way to keep frequently used supplies at your fingertips. Stable base allow for easy positioning for work or for moving out of the way when not in use.

Phlebotomy Stands & Work Stations feature:

- Cleanable work surfaces
- Adjustable heights
- Durable construction



Two-Bin, Mobile, Phlebotomy Stand

- ClintonClean™ solid plastic top
- Chrome plated steel pole with knob-lock adjustment
- 5-leg, cast aluminum base with bumper guards
- 2" dual wheel, soft roll casters
- Two easy to remove drop-in bins/11" L x 8¹/₄" W x 3¹/₂" D
(28 cm L x 21 cm W x 9 cm D)

6902

top	width 13 ¹ / ₂ " 34.29 cm	depth 21" 53.34 cm
base diameter	23" 58.42 cm	
height range	30" – 38" 76.2 cm – 96.52 cm	

(Supplies shown not included)



6902

One-Bin, Mobile, Phlebotomy Stand

- ClintonClean™ solid plastic top
- Chrome plated, steel pole with knob-lock adjustment
- 4-leg, epoxy-coated stand with bumper guards
- 2" rubber wheel, ball bearing casters
- One easy to remove drop-in bin/16¹/₂" L x 11" W x 2¹/₂" D
(42 cm L x 28 cm W x 6 cm D)

6901

top	width 15" 38 cm	depth 21 ¹ / ₂ " 54.6 cm
base diameter	22" 55.88 cm	
height range	30" – 38" 76.2 cm – 96.52 cm	



6901

Mobile, Phlebotomy Work Station

- Easy-clean, vacuum-formed, gray, ABS, plastic top with textured finish and raised rim on all sides
- 5-leg, cast aluminum base with bumper guards
- 2" dual wheel, soft roll casters

6900

top	width 15 ³ / ₄ " 40 cm	depth 18" 45.72 cm
base diameter	24" 60.96 cm	
height range	29 ¹ / ₂ " – 36" 74.93 cm – 91.44 cm	



6900

(Supplies shown on all models not included)

Mobile, Phlebotomy Work Station with Bin

- Easy-clean, vacuum-formed, gray, ABS, plastic top with textured finish and raised rim on all sides
- 5-leg, cast aluminum base with bumper guards
- 2" dual wheel, soft roll casters
- One easy to remove drop-in bin/10³/₄" L x 8¹/₄" W x 3¹/₂" D
(28 cm L x 21 cm W x 9 cm D)

6900-B

top	width 27 ¹ / ₄ " 70.48 cm	depth 18" 45.72 cm
base diameter	24" 60.96 cm	
height range	27 ¹ / ₄ " – 42 ¹ / ₄ " 69.21 cm – 107.31 cm	



6900-B

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.