

Mobile Phlebotomy Work Station with Bin

MODEL 6900-B

- Vacuum-formed, ABS plastic top with textured finish
- Top is rimmed on all 4 sides
- 5-leg, cast aluminum base
- 2" dual wheel, soft roll casters
- One 10-3/4" L x 8-1/4" W x 3-1/2" D Drop in plastic bin (27cm L x 21cm W x 9cm D)

NOTE: *Supplies pictured not included.*

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.



MODEL 6900-B SPECIFICATIONS

Top 27.25"L x 18"W

Base 24" D

Height Range 27.25"— 42.25"

Prop 65 Warning



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1